MeRA Test

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MeRA Test

ENGLISH

Microbiological test with spores of *Geobacillus stearothermophilus*
for the detection of antibiotic and sulphamide residues in meat

DESCRIPTION

Antimicrobial substances are given to cattle for therapeutic treatment of infections and so they can be present in meat as residual drugs. The presence of residual antimicrobial drugs in meat is a potential hazard for the consumers since those substances interfere with the intestinal flora, can cause allergic reactions and can contribute to the spreading of resistant microorganisms.

The detection of inhibiting substances in meat is ruled in Europe by Regulation 37/2010 EC and following modifications.

MeRA (Meat Rapid Antibiotic) Test is a microbiological test, set for the detection of residual antibiotics and sulphamides in meat, that uses spores of *Geobacillus stearothermophilus* and a culture medium designed for that purpose.

CONTENT OF THE PACKAGE

Each package contains:

- 50 tubes containing liquid culture medium
- 50 discs containing spores of *Geobacillus stearothermophilus*
- 1 instructions sheet

PRINCIPLE OF THE METHOD

The test detects residues of antibiotics and sulphamides through their inhibition action against the growth *in vitro*, at $64 \pm 0,5$ °C, of *Geobacillus stearothermophilus*. The inhibition of bacterial growth is made evident by the presence of a fermentable sugar and a pH indicator in the culture medium.

PROCEDURE FOR USE

- Weigh about g 50 of meat. Cut up the tissue and mince further by electric mixer.
- Weigh g 2 of minced meat and introduce it in a 10 mL tube.
- Add mL 6 of distilled water to the tube (ratio meat:water 1:3).
- Homogenize for some seconds through homogenizer set at 13000 rpm.
- Centrifuge the homogenized meat by tabletop centrifuge set at 4000 rpm for 15 minutes.
- Collect the supernatant (meat extract) and make sure its pH is 5.6 roughly. Otherwise, add NaOH 0.1 M or HCl 0.1 M respectively to alkalize or acidify the meat extract.

NOTE: pH lower than 5.6 causes the immediate color change to yellow of MeRA Test thus invalidating the test; pH above 5.6 causes delay of the color change in case of antibiotic free meat sample.

- Remove the cup from the tube and add 1 disc of spores. Shake gently.
- Put the cup and introduce the tube in a waterbath at $64 \pm 0,5$ °C or in a Termoblock set at $64 \pm 0,5$ °C. Incubate for 20'.
- Remove the tube from the waterbath or Termoblock and allow to reach room temperature (roughly 5').
- Remove the cup and introduce 1 mL of the supernatant under test.
- Prepare one or more analyte-positive controls, by the addition of antibiotic/sulphamide into the tubes as indicated in the paragraph "Quality Control for the user".
- Replace and tighten the cap on the tubes and let the antibiotic/sulphamide, if present, act at room temperature for 20'.
- Reintroduce the tube in the waterbath or in the Termoblock at $64 \pm 0,5$ °C for the second incubation for 3h - 3h 30'.

INTERPRETATION OF THE RESULTS

At the end of the second incubation observe colour of the liquid in the tubes that, if green/blu, indicates antibiotic/sulphamide residues in the examined meat in concentrations above the detection limits; otherwise, a yellow/orange colour indicates the absence of antibiotic/sulphamide residues in the examined meat or their presence in concentrations under the detection limits.

QUALITY CONTROL FOR THE USER

Every batch of **MeRA Test** is submitted to quality control. The analyte-negative control is obtained using a certified antibiotic free meat sample. The analyte-positive controls are obtained using meat samples with the addition of antibiotics/sulphamides, such as ampicillin (50 µg/kg) and sulphadiazine (100 µg/kg).

PRECAUTIONS

The product **MeRA Test** is not classified as dangerous under current legislation. **MeRA Test** is a disposable device to be utilized only for *in vitro* use, it is intended for use in a professional environment and must be used in laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store **MeRA Test** at 2-8 °C in its original packaging. Keep away from heat sources and avoid excessive temperature changes. In such conditions **MeRA Test** will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using it if there are deterioration signs.

DISPOSAL

After use, before disposal, tubes must be properly decontaminated with chemical substances or autoclaved.

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Table 1: Sensitivity of **MeRA Test** to some antimicrobial agents, commonly used in veterinary medicine and respective MRL (Maximum Residue Limits) values valid for Europe¹.

Table 1

ANTIBIOTIC/SULPHAMIDE	MRL for bovine muscle ¹ MAXIMUM RESIDUE LIMITS (µg/kg)	MeRA Test Sensitivity in 3h 30' DETECTION LIMITS (µg/kg)
BETA-LACTAMS		
Ampicillin	50	25-50
Oxacillina	300	150-300
Cloxacillin	300	150-300
Dicloxacillina	300	150-300
Amoxycillin	50	25-50
Benzylpenicillin	50	25-50
Penethamate	50	25-50
Cefalexin	200	100-200
Ceftiofur	1000	<500
Cefquinome	50	25-50
TETRACYCLINES		
Tetracycline	100	100-200
Chlortetracycline	100	100-200
Oxitetracycline	100	100-200
Doxycycline	100	100-200
MACROLIDES		
Erythromycin	200	200-400
Tilosina	100	100-200
Tilmicosin	50	100-200
Spiramicina	200	200-400
LINCOSAMIDES		
Lincomycin	100	100-200
Pirlimycin	100	50-100
AMINOGLYCOSIDES		
Gentamicin	50	100-200
Neomycin	500	500-1000
Streptomycin	500	500-1000
Dihydrostreptomycin	500	500-1000
SULPHAMIDES		
Sulfadiazine	100	50-100
SULPHANILAMIDES		
Sulfadimidine	100	50-100
BENZYL PYRIMIDINE		
Trimethoprim	50	50-100
QUINOLONES		
Flumequine	200	200-400
Enrofloxacin	100	50-100

¹ Regulation 37/2010 EC.

TABLE OF SYMBOLS

 Batch code	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	

Liofilchem®

Via Scozia Zona Ind.le - 64026 Roseto D.A. (TE) - Italy

Tel. +390858930745 Fax +390858930330 Website: www.liofilchem.net E-Mail: liofilchem@liofilchem.net

F00041
Rev.5.1 / 03.02.2016