
# REF M42 MICROGEN™ SALMONELLA

## INTENDED USE

Microgen™ Salmonella is a latex slide agglutination test for the confirmatory identification of presumptive Salmonella colonies from selective agar plates. The kit is intended for professional use only.

## PRINCIPLE OF THE TEST

Latex particles are coated with polyvalent antisera raised against a wide range of Salmonella antigens. When mixed with a suspension of Salmonella organisms, the latex particles rapidly agglutinate to form visible clumps. Microgen™ Salmonella detects >99% of motile Salmonella species and early investigations have indicated that specific non-motile species may also be detected.

CONT		KIT PRESENTATION	
REAG	TEST	M42a	Salmonella Latex Reagent: 2.5mL
Latex particles coated with rabbit antiserum against Salmonella antigens. Preserved with 0.099% sodium azide. (Blue cap)			
CONTROL	+	M42b	Positive Control: 0.5mL
Inactivated preparation of Salmonella antigens preserved with 0.099% sodium azide. (Black cap)			
NaCl	0.85%	M40	0.85% Isotonic Saline: 5.0mL
Preserved with 0.099% sodium azide. (White cap)			

Instructions for Use  
Disposable agglutination slides  
Disposable mixing sticks

## Additional Requirements:

- Bacteriological loops

## WARNINGS AND PRECAUTIONS

### Safety:

- The reagents supplied in this kit are for *in vitro* diagnostic use only
- Sodium azide, which is used as a preservative in the kit reagents can react with lead or copper plumbing to form potentially explosive metal azides. Dispose by flushing with a large volume of water to prevent azide build-up.
- Appropriate precautions should be taken when handling or disposing of potential pathogens. Decontamination of infectious material can be achieved with sodium hypochlorite at a final concentration of 3% for 30minutes. Liquid waste containing acid must be neutralised before treatment.
- The positive control has been inactivated during the manufacturing process. However, it should be handled as though potentially infectious.

## Procedural:

- Microgen™ Salmonella should be used according to the kit instructions.
- Allow all reagents to reach room temperature before use.
- Do not dilute any of the kit reagents
- Do not intermix reagents from different batches of kits.
- Do not freeze any of the kit reagents
- Do not allow the latex reagent dropper to touch the positive control or bacterial samples.
- Be careful only to record agglutination. Reactions that are “curdy” or “stringy” may not be true agglutination.
- Ensure the slide is clean and dry prior to use.

## STORAGE AND SHELF LIFE

Microgen™ Salmonella should be stored at 2-8°C when not in use. The kit should not be used after the expiry date printed on the carton label.

## SPECIMENS

Colonies grown on selective agar plates can be tested with Microgen™ Salmonella.

## PROCEDURE

### Quality Control:

The following controls should be performed each time the kit is used.

- Reagent Control: Add one drop of Microgen™ Salmonella latex (M42a) to one drop of M40 saline solution in the same circle on a slide. Mix and spread the liquid over the entire area of the circle with a mixing stick. Rock the slide gently for 2 minutes and observe for agglutination. If any agglutination is seen, either the latex or the saline is contaminated and should be discarded.
- Positive Control: Add one drop of positive control (M42b) to one circle on the test slide. Add one drop of Microgen™ Salmonella latex to the same circle and mix. **Do not allow the dropper to touch the positive control.** Rock the slide gently. Within 2 minutes, agglutination, indicating a positive result, should be visible. If no agglutination is seen, a fresh kit should be used.

### Test Procedure:

- Dispense 1 drop of M40 isotonic saline into a circle of a Microgen™ agglutination slide.
- Using an inoculating loop, remove a colony from the selective agar plate and emulsify the colony in the drop of saline to produce a heavy smooth suspension. Suspensions should only be made from colonies with morphologies resembling *Salmonella spp.*
- Rock the slide gently for up to 2 minutes and observe for autoagglutination or clumping. If the suspension remains smooth, proceed to Step 4 (see Limitations of Use Note 1).
- Mix the Microgen™ Salmonella latex by gently inverting and add one drop next to the bacterial suspension. **Do not allow the dropper to touch the suspension.**

- Mix the latex reagent and the bacterial suspension with a clean mixing stick and rock the slide gently two or three times. Excessive rocking of the slide is not necessary. Examine for agglutination within a maximum of 2 minutes.
- After reading, discard the used mixing sticks and slides into suitable disinfectant.

## INTERPRETATION

Agglutination within 2 minutes is a positive result and indicates the presence of *Salmonella* in the sample. Absence of agglutination indicates that *Salmonella* is not present in the test culture.

## LIMITATIONS OF USE

- Results should be interpreted in the context of all available clinical and laboratory information.
- Rough strains of *Salmonella* are known to cause non-specific autoagglutination in saline alone and therefore cannot be tested with Microgen™ *Salmonella*.
- Some non-motile strains may not be detected by Microgen™ *Salmonella*.
- Some oxidase-positive organisms may give false positive reactions.
- Old stock cultures of Enterobacteriaceae on nutrient agar slopes may cause non-specific agglutination whereas old stocks of *Salmonella* may give false negative results. Fresh sub-cultures should be prepared for testing.
- Identification with Microgen™ *Salmonella* is presumptive and all positive results should be confirmed by further identification tests and serotyping of pure cultures.

## PERFORMANCE CHARACTERISTICS

Microgen™ *Salmonella* has been evaluated in comparison with a well-established commercially available latex agglutination test for *Salmonella*. 126 isolates of *Salmonella* spp. and a range of 58 potentially cross-reacting bacteria were tested in both products.

		Microgen™ <i>Salmonella</i>		Total
		+ve	-ve	
Commercial Latex Test	+ve	135**	1*	136
	-ve	0	48***	48
Total		135	49	184

Sensitivity:  $135/136 = 99.3\%$

Specificity:  $48/48 = 100\%$

Concordance:  $183/184 = 99.5\%$

\*1 sample was negative in Microgen™ *Salmonella* but equivocal in the commercial test. This sample was subsequently identified as *Salmonella bergien*.

\*\* Of the 135 isolates in this group, 11 were cross reactants in both tests. These were isolates of *C. diversus* (1), *A. baimannii* (2), *P. stuartii* (1), *B. cereus* (2), *S. aureus* (4), *Strep spp* (1)

However all of the above either did not grow or showed very atypical morphologies, when cultured on *Salmonella*-selective media. In the case of *B. cereus*, agglutination was atypical (stringy)

\*\*\* 1 *S. dublin* was repeatedly negative in both tests.

## REPRODUCIBILITY

**Intra-batch reproducibility** was evaluated by testing sensitivity and specificity of 1 batch of product against serial dilutions of reference and kit control antigens, and a panel of 34 bacterial samples. Different operators carried out tests on 3 separate occasions. End-point titres obtained with reference/control antigens and qualitative results with the panel were identical in the three assays.

**Inter-batch reproducibility** was examined by testing sensitivity and specificity of 3 batches of product against serial dilutions of reference and kit control antigens, and a panel of 34 bacterial samples. Between the 3 batches, variation in end-point titres was minimal (1

doubling dilution) and qualitative results with the panel correlated 100%.

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Microgen Bioproducts Ltd