
REF M45 MICROGEN® LEGIONELLA

INTENDED USE

Microgen® Legionella is a latex agglutination test intended for confirmatory identification of *Legionella pneumophila* and commonly isolated *Legionella* species grown on selective media. The test is suitable for organisms derived from patients with suspected Legionella pneumonia or from environmental sources. Microgen® Legionella allows the separate identification of *L. pneumophila* Serogroup 1 and Serogroups 2-15* and commonly isolated *Legionella* species.

*Note: *L. pneumophila* serogroup 16 is now mentioned in the Manual of Clinical Microbiology. Serogroup 16 should cross-react with serogroup 6, so theoretically the Microgen 2-15 latex reagent should also agglutinate with serogroup 16 isolates.

The kit is intended for professional use only.

PRINCIPLE OF THE TEST

For Test Reagent 1, latex particles are coated with polyclonal rabbit antibodies raised against *L. pneumophila* Serogroup 1. For Test Reagent 2-15, latex particles are coated with polyclonal rabbit antibodies raised against serogroups 2-15. For Test Reagent Species, latex particles are coated with polyclonal rabbit antibodies raised against 10 commonly isolated *Legionella* species (see table below). When these latex particles are mixed with a suspension containing the appropriate *Legionella* bacteria or heat killed antigens from the relevant *Legionella* bacteria, an immunochemical reaction takes place causing the latex particles to agglutinate into aggregates which are easily visible to the naked eye. The latex reagent coated with serogroup 1 antibodies will only agglutinate in the presence of serogroup 1 antigens. The latex reagent coated with serogroup 2-15 antibodies will only agglutinate in the presence of antigens from any one of the serogroups 2-15. The latex reagent coated with *Legionella* species antibodies will agglutinate with the 10 relevant *Legionella* species (see table below)

Reactions of Test Reagent 1 with serogroup 1 antigens are generally stronger and faster than with those between Test Reagent 2-15 and serogroup 2-15 antigens. This is because of the dilution effect of blending 15 different antisera during manufacture of Test Reagent 2-15. The *Legionella* species latex reacts in a similar manner to the *Legionella pneumophila* serogroup 1 latex.

CONT

KIT PRESENTATION

REAG	TEST	1	M45a Test Reagent 1:	2.5mL
------	------	---	----------------------	-------

Latex particles coated with rabbit polyclonal antibodies to *L. pneumophila* serogroup 1.
Preserved with 0.099% sodium azide. (Blue cap)

REAG	TEST	2-15	M45b Test Reagent 2-15:	2.5mL
------	------	------	-------------------------	-------

Latex particles coated with polyclonal rabbit antibodies to *L. pneumophila* serogroups 2-15. Preserved with 0.099% sodium azide. (Yellow cap)

CONTROL	+
---------	---

M45c Positive Control: 1.0mL

Suspension of inactivated *Legionella* antigens reactive with Test Reagents 1, 2-15 and Species. Preserved with 0.099% sodium azide. (Black cap)

REAG	TEST	SPECIES
------	------	---------

M45d Test Reagent Species 2.5mL

Latex particles coated with polyclonal rabbit antibodies to 10 commonly isolated *Legionella* species (see table below)
Preserved with 0.099% sodium azide (Red cap)

NaCl	0.85%
------	-------

M40 0.85% Isotonic Saline: 2 x 5.0mL

Preserved with 0.099% sodium azide. (White cap)

Instructions for Use
Disposable agglutination slides
Disposable mixing sticks

Additional Requirements:

- Bacteriological loops
- Legionella selective medium
- Glass tubes for boiling
- 0.85% isotonic saline

WARNINGS AND PRECAUTIONS

Safety:

- The reagents supplied in this kit are for *in vitro* diagnostic use only
- Sodium azide, which is used as a preservative in the kit reagents, can react with lead or copper plumbing to form potentially explosive metal azides. Dispose by flushing with a large volume of water to prevent azide build-up.
- Appropriate precautions should be taken when handling or disposing of potential pathogens. Decontamination of infectious material can be achieved with sodium hypochlorite at a final concentration of 3% for 30 minutes. Liquid waste containing acid must be neutralised before treatment.
- The positive control has been inactivated during the manufacturing process. However, it should be handled as though potentially infectious.

Procedural:

- Microgen® Legionella should be used according to the kit instructions.
- Allow all reagents to reach room temperature before use.
- Do not dilute any of the kit reagents
- Do not intermix reagents from different batches of kits.
- Do not freeze any of the kit reagents
- Do not allow the latex reagent dropper to touch the positive control or bacterial samples.
- Be careful only to record agglutination. Reactions that are "curdy" or "stringy" may not be true agglutination.
- Ensure the agglutination slide is clean and dry prior to use.

STORAGE AND SHELF LIFE

Microgen® Legionella should be stored at 2-8°C when not in use. The kit should not be used after the expiry date printed on the carton label.

SPECIMENS

Colonies grown on selective agar plates can be tested with Microgen® Legionella. The morphology of the colonies tested should resemble that of *Legionella*.

PROCEDURE

Quality Control:

The following controls should be performed each time the kit is used.

1. Reagent Control: Gently mix the Microgen® test latex reagents (M45a, M45b, M45d) and add one drop of each reagent to separate wells on the agglutination slide. Add one drop of saline solution (M40) to each drop of latex reagent. Using a different mixing stick for each well, mix the latex reagent and saline thoroughly, spreading the liquid over the entire area of the well. Rock the slide gently for 2 minutes and observe for agglutination. If any agglutination is seen, at least one of the reagents is contaminated and a fresh kit should be used.
 2. Positive Control: Add one drop of positive control (M45c) to one well on the agglutination slide. Gently mix Test Reagent 1 (M45a) and add 1 drop to the same well. **Do not allow the dropper to touch the positive control.** Mix with a mixing stick, rock the slide gently for 2 minutes and observe for agglutination. A positive result, indicated by agglutination, should be seen. Repeat this process using Test Reagent 2-15 (M45b) and Test Reagent Species (M45d). Again, a positive result should be seen. If agglutination is not seen with any of the three latex reagents, a fresh kit should be used.
- NOTE: The reaction strength with the positive control may not be the same for Test Reagent 1 as with Test Reagent 2-15 and Test Reagent Species. (see PRINCIPLE OF THE TEST above)

Test Procedure:

1. Dispense 1 drop of isotonic saline (M40) on to each of three wells of the agglutination slide.
2. Using an inoculating loop, remove 3-4 Legionella-like colonies from the selective agar plate and make a thick even smear on the slide alongside each drop of saline.
3. Mix the colonies with the saline and emulsify to form a smooth heavy suspension, spreading the liquid over the entire surface of the well.
4. If the suspension remains smooth, proceed to section 7. If the suspension is "stringy" or "granular" (often as a result of old and/or mucoid cultures), proceed as follows:
5. Dispense 0.5mL 0.85% isotonic saline into a glass tube. Prepare a homogeneous turbid suspension of organisms taken from the selective agar plate.
6. Boil the suspension for 5 minutes. Allow to cool to room temperature. Place 30µL boiled suspension on to each of three wells of an agglutination slide.
7. Gently mix each test latex reagent (M45a, M45b and M45d) to ensure a homogeneous suspension
8. Add one drop of Test Reagent 1 to one of the bacterial suspensions, one drop of Test Reagent 2-15 to the second suspension and one drop of Test Reagent Species to the third suspension. **Do not allow the reagent dropper to touch the suspension.**
9. Mix the reagent and suspension using a new mixing stick for each combination. Spread the liquid over the entire area of the well.
10. Rock the slide gently for 2 minutes and observe for agglutination.
11. Discard the used mixing sticks and slides into a suitable disinfectant.

INTERPRETATION

An agglutination reaction is indicated by visible aggregation of the latex particles. Microgen® Legionella results should be interpreted as follows:

Test Reagent 1	Test Reagent 2-15	Test Reagent Species	Interpretation
+	-	-	<i>L. pneumophila</i> 1 present
-	+	-	<i>L. pneumophila</i> 2-15 present
-	-	+	<i>Legionella</i> species present
-	-	-	No <i>Legionella</i> present
+	+	+	Possible non specific agglutination, inconclusive result*
+	+	-	Possible non specific agglutination, inconclusive result*
+	-	+	Possible non specific agglutination, inconclusive result*
-	+	+	Possible non-specific agglutination. Inconclusive result*

*A non-specific agglutination pattern does not preclude the presence of *Legionella* bacteria but results have to be interpreted as inconclusive.

LIMITATIONS OF USE

1. Results should be interpreted by the clinician in the context of all available clinical and laboratory information.
2. Stringy reactions on the slide may not be true positive reactions and further tests are required.
3. Old and/or mucoid cultures may not give a smooth suspension in saline and may give atypical agglutination. These should be pre-treated by boiling as described in the method above.
4. Microgen® Legionella is intended for the identification of *L. pneumophila* serotypes 1-15 and detection of a range of 10 *Legionella* species following culture on selective agar plates. Colonies giving positive results should be confirmed as *Legionella* spp. by their inability to grow on Legionella selective culture media deficient in cysteine.

PERFORMANCE CHARACTERISTICS

Microgen® Legionella has been evaluated in comparison with a well-established, commercially available latex agglutination test for *L. pneumophila*. 128 isolates of *L. pneumophila*, non-pneumophila *Legionella* spp and potentially cross-reacting bacteria were tested in both products.

		Microgen® Legionella (all serogroups)		Total
		+ve	-ve	
Commercial Latex Test (all serogroups)	+ve	50*	0	50
	-ve	0	78	78
Total		50	78	128

Sensitivity: 50/50 = 100%

Specificity: 78/78 = 100%

Accuracy: 128/128 = 100%

* Of the 50 isolates in this group, 12 were cross reactants in both tests. These were isolates of *S. aureus* (4), *C. diversus* (1), *A. baumannii* (2), *P. stuartii* (1), *B. cereus* (2), *K. pneumoniae* (1), *Streptococcus* spp. (1). However all of the above either do not grow or show very atypical morphology, when cultured on Legionella selective media. In the case of *B. cereus*, agglutination was atypical (stringy)

70 isolates of *L. pneumophila* and *Legionella spp* were tested in both Microgen® Legionella and the comparison latex kit to differentiate between Serogroup 1 and Serogroups 2-15.

		Microgen® Legionella			
		SG1 Positive	SG2-15 Positive	L. pn. Negative	Total
Commercial Latex Test	SG1 Positive	8	0	0	8
	SG2-15 Positive	0	30	0	30
	L. pn. negative	0	0	32	32
	Total	8	30	32	70

No isolates of *L. pneumophila* serogroup 1 cross-reacted with Test Reagent 2-15 and, similarly, no isolates of *L. pneumophila* serogroups 2-15 cross-reacted with Test Reagent 1. All non-pneumophila *Legionella spp* isolates were non-reactive with Test Reagent 1 and Test Reagent 2-15 in Microgen Legionella.

Legionella Species Latex has been tested against all *Legionella pneumophila* serogroups with no cross-reactivity detected. The Legionella Species Latex is designed to detect 10 *Legionella* species. In particular these are:

Legionella species (10 Target Species)
<i>L. micdadei</i>
<i>L. bozemani</i> 1
<i>L. bozemani</i> 2
<i>L. dumoffi</i>
<i>L. longbeachae</i> 1
<i>L. longbeachae</i> 2
<i>L. jordanis</i>
<i>L. gormanii</i>
<i>L. anisa</i>
<i>L. feelei</i>

Specificity testing of the Legionella Species Latex indicates that isolates of *Legionella sainthelensi* sgp 1, *erythra*, *hacklaie* sgp 1 and *tucsonensis* can also produce positive results.

REPRODUCIBILITY

Intra-batch reproducibility was established by testing sensitivity and specificity of 1 batch of product against serial dilutions of reference and kit control antigens and a panel of 47 bacterial samples. Different operators carried out tests on 3 separate occasions. End-point titres obtained with reference/control antigens and qualitative results with the panel were identical in the three assays.

Inter-batch reproducibility was examined by testing sensitivity and specificity of 3 batches of product against serial dilutions of reference and kit control antigens, and a panel of 47 bacterial samples. Between the 3 batches, no significant differences in end-point titres were seen and qualitative results with the panel correlated 100%.

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Microgen Bioproducts Ltd