## THAYER MARTIN AGAR

Selective medium for *Neisseria spp* isolation.

### TYPICAL FORMULA (g/L)

Proteose Peptone N°3.....	15.0
Corn Starch.....	1.0
Potassium Phosphate dibasic.....	4.0
Potassium Phosphate monobasic.....	1.0
Sodium Chloride.....	5.0
Agar.....	17.0
Horse defibrinated blood.....	50.0 ml
Cyanocobalamine.....	0.05 mg
Glutamine.....	5.5 mg
Thiamine HCl.....	1.0 ml
L-Cysteine HCl.....	0.12 ml
Adenine.....	5.0 mg
p-aminobenzoic acid.....	0.065 mg
Coccarboxylase.....	0.5 mg
Guanine.....	0.15 mg
NAD.....	1.25 mg
Ferric nitrate.....	0.1 mg
Vancomycin.....	1.0 mg
Colistin sulphate.....	3.75 mg
Amphotericin B.....	0.5 mg
Trimethoprim.....	1.5 mg
Final pH=7.2 ± 0.2	

### DESCRIPTION

THAYER MARTIN AGAR is used to isolate and cultivate *Neisseria spp*.

### PRINCIPLE

*Neisseria gonorrhoeae* and *Neisseria meningitidis* are aerobic, Gram-negative diplococci, which grow optimally in an environment with 2-10% CO<sub>2</sub>. These organisms are very fastidious and cultivate with oxidase positive colonies.

### TECHNIQUE

Inoculate by rolling a culture swab across a segment of the plate or in a large "Z" pattern so that an adequate area of the plate is inoculated. Streaking of the plate is carried out with a sterile loop to ensure an adequate dispersion of the organisms. Plates are incubated at 36±1 °C, in a sealed jar with at least 70% humidity and 5-10% carbon dioxide.

Examine the plates after 24 hours incubation and, if negative, reincubate for further 24 hours.

Presumptive gonococcus colonies are identified by the Gram stain, oxidase and sugar fermentation reactions.

### INTERPRETATION of RESULTS

Typical colonies isolated on this medium which are oxidase positive and composed of Gram-negative diplococci can be considered to be presumptive of *Neisseria gonorrhoeae* and *Neisseria meningitidis*.

### STORAGE

2-8 °C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

### WARNING and PRECAUTIONS

For laboratory use.

The product doesn't contain dangerous substances according to directives 1999/45/CE and 2001/60/CE or for which exist recognized exposure limits.

### DISPOSAL of WASTE

Disposal of waste must be carried out according to national and local regulations in force.

### REFERENCES

1. Chapin, C.K., G.V. Doern (1983). *J. Clin. Microbiol.* 17: 1163-1165.
2. Martin, J.E., Armstrong J.H., Smith P.B. (1974). *Appl. Microbiol.* 27: 802-805.
3. NCCLS document M22-A2, 1996. *Quality Assurance for Commercially prepared Microbiological Culture Media-Second ed. Approved Standard.*
4. Seth, A. (1970). *Brit. J. Vener. Dis.* 46: 201-202.

**Liofilchem s.r.l** Bacteriology Products

Via Scozia-Zona industriale - 64026 Roseto degli Abruzzi Tel. +39.085.8930745 - Fax +39.085.8930330  
Web site: <http://www.liofilchem.net> E-mail: [liofilchem@liofilchem.net](mailto:liofilchem@liofilchem.net)

## PRODUCT SPECIFICATIONS

### NAME

THAYER MARTIN AGAR

### PRESENTATION

Ready-to-use plates containing 22+/-1 ml.

### PACKING

Code	Content	Packaging
11040	20 plates	<ul style="list-style-type: none"> <li>5 plates in thermically soldered film</li> <li>4 x 5 plates in cardboard boxes</li> </ul>
11040*	100 plates	<ul style="list-style-type: none"> <li>5 plates in thermically soldered film</li> <li>2 x 5 plates in thermically soldered bag</li> <li>10 piles (2 x 5 ps) in box</li> </ul>

### pH OF THE MEDIUM

7.2 ± 0.2.

### STORAGE

2-8 °C.

### USE

THAYER MARTIN AGAR is used to isolate and cultivate *Neisseria spp.*

### TECHNIQUE

Refer to technical sheet of the product.

### APPEARANCE OF THE MEDIUM

Opaque chocolate-brown medium.

### SHELF LIFE

4 months

### QUALITY CONTROL

- Control of general characteristics, label and print.
- Sterility control:  
7 days at 25 ± 1 °C, in aerobiosis.  
7 days at 36 ± 1 °C, in aerobiosis.
- Microbiological control:  
Inoculum for productivity: 10-100 UFC/ml.  
Inoculum for selectivity: 10<sup>4</sup>-10<sup>5</sup> UFC/ml.  
Inoculum for specificity: ≤ 10<sup>4</sup> UFC/ml.  
Incubation conditions: 24-48 h at 36 ± 1 °C / 5-10% CO<sub>2</sub>.

Microorganism		Growth
<i>Neisseria gonorrhoeae</i>	ATCC 43069	Good
<i>Neisseria meningitidis</i>	ATCC 13090	Good

### TABLE of SYMBOLS

Symbol	Meanings
	Catalogue number
	In vitro Diagnostic Medical Device
	Manufacturer
	Temperature limitation
	Kit content
	Use by
	Batch code
	Consult accompanying documents

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