Hektoen Enteric Agar

Selective and differential medium for detection of pathogenic intestinal bacteria from clinical and nonclinical specimens.

INTENDED PURPOSE

Selective medium for the isolation and presumptive identification of *Salmonella* and *Shigella* spp. from clinical and nonclinical specimens. This medium is intended as an aid in the diagnosis, requiring further identification tests to complete the diagnostic results.

DESCRIPTION

Hektoen Enteric (HE) Agar is a moderately selective medium used for the isolation and cultivation of Gram-negative enteric microorganisms, especially *Salmonella* and *Shigella*, from faeces, rectal swabs, food, animal feed and environmental materials.

This medium complies to ISO 21567 and FDA-BAM for detection of *Shigella* and *Salmonella*, respectively.

TYPICAL FORMULA*

	(g/litre)
Enzymatic Digest of Meat	12.0
Yeast Extract	3.0
Lactose	12.0
Saccharose	12.0
Salicin	2.0
Bile Salts No. 3	9.0
Sodium Chloride	5.0
Sodium Thiosulfate	5.0
Ammonium Ferric Citrate	1.5
Acid Fuchsin	0.1
Bromthymol Blue	0.065
Agar	15.0
Final pH 7.5 ± 0.2 at 25°C	

*Adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

Enzymatic digest of meat provides amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of B-group. Lactose, saccharose and salicin are fermentable carbohydrates. Bile salts and acid fuchsin inhibit Gram-positive organisms. Sodium chloride maintains the osmotic balance of the medium. Ammonium ferric citrate and sodium thiosulfate enable the detection of hydrogen sulfide production. Bromothymol blue together with acid fuchsin act are the pH indicator system. Agar is the solidifying agent.

PREPARATION

Dehydrated medium

Suspend 76 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. DO NOT AUTOCLAVE.

Medium in bottles

Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: Autoclave, water bath, sterile Petri plates, test tubes, inoculating loops, swabs, incubator, quality control organisms.

SPECIMENS

Faecal specimens should be sampled at the acute stage, before antimicrobial therapy (where possible) and examined as soon as possible after collection. Ensure that specimen is placed in appropriate package (clean,

leak-proof container) and, in case the laboratory processing is delayed (>1 h), use a suitable transport medium, such as Cary-Blair for transportation/storage. Refer to specific guidelines for more information about specimen collection and preparation.

TEST PROCEDURE

Inoculate the plates by directly streaking the specimen on the agar surface or spread the sample from an enrichment culture to obtain well-isolated colonies.

Incubate aerobically at $35 \pm 2^\circ\text{C}$ for 18-24 h.

For more details, consult appropriate guidances.

INTERPRETING RESULTS

Shigella and *Providencia* spp, form green, moist colonies.

Salmonella and *Proteus* spp, grow as blue-green colonies, with or without black center due to H_2S production.

Coliforms, which are mostly rapid lactose-saccharose-salicin fermenters, develop red-salmon colonies surrounded by a zone of bile precipitate.

Enterococci, Staphylococci and other Gram-positive bacteria are partially or completely inhibited.

Note: further testing should be conducted to confirm the presumptive identification of organisms isolated on this medium.

STORAGE

The powder is very hygroscopic, store the powder at $10-30^\circ\text{C}$, in a dry environment, in its original container tightly closed. Store bottles and prepared plates at $10-25^\circ\text{C}$ away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.

Medium in bottles: 2 years.

Ready-to-use plates: 6 months.

QUALITY CONTROL

Appearance of Dehydrated Medium: Free-flowing, homogeneous, light beige.

Appearance of Prepared Medium: Slightly opalescent, greenish to brownish.

Expected Cultural Response:

Control strain	Inoculum	Incubation	Criteria	Specification
<i>Salmonella</i> Typhimurium	WDCM 00031 (ATCC® 14028; NCTC 12023)	50-100 CFU	Good growth ($P_R \geq 0.5$)	Blue-green colonies with black centre
<i>Shigella flexneri</i>	WDCM 00126 (ATCC® 12022; NCTC 12698)			Green colonies
<i>Proteus mirabilis</i>	ATCC® 12453	10^3-10^4 CFU	Weak to good growth	Blue-green colonies with black centre
<i>Klebsiella pneumoniae</i>	WDCM 00097 (ATCC® 13883; NCTC 9633)			Red-salmon colonies with zone of bile precipitate
<i>Escherichia coli</i>	WDCM 00012 (ATCC® 8739; NCTC 12923)			Red-salmon colonies with or without zone of bile precipitate
<i>Enterococcus faecalis</i>	WDCM 00087 (ATCC® 29212; NCTC 12697)	10^4-10^6 CFU	None or poor growth	—

A productivity ratio (P_R) of 0.5 is equivalent to a recovery rate of 50%.

Please refer to the actual batch related Certificate of Analysis (CoA).

PERFORMANCE CHARACTERISTICS

Performance testing of HE agar was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Invalid results can be caused by poor specimen quality, improper sample collection, improper transportation, improper laboratory processing, or a limitation of the testing technology. The operator should understand the principles of the procedures, including its performance limitations, in advance of operation to avoid potential mistakes.

Due to nutritional variation, some strains may result in poor growth or fail to grow on this medium. Additional agar media, such as MacConkey agar, should be used alongside the HE agar for primary isolation.

Colonies of *Proteus* may resemble *Salmonella* or *Shigella*. Rare strains of *Salmonella* and *Shigella* are capable of lactose fermentation, producing red-salmon colonies which may resemble coliforms.

HE agar is intended as an aid in the diagnosis of infectious diseases, requiring further tests to complete the diagnostic results. All identification tests should ideally be performed from non-selective agar.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

The product is available in the various configurations listed below. There may be additional product ref. numbers as well. For an updated listing of available products, visit liofilchem.com


Product	Format	Packaging	Ref.
Hektoen Enteric Agar	Plate 90 mm	20 plates	10043
		100 plates	10043*
	Bottle	6 x 100 ml	402230
		6 x 200 ml	412230
	Dehydrated media	100 g	620021
		500 g	610021
		5 kg	6100215

Significant changes from previous version (rev.0, 2015-08-04):

Document	Release Date	Change Summary
10043_IFU-1	2022-05-16	Updated: Layout and content in compliance with IVDR 2017/746 Added: Example photographs

This IFU document and the SDS are available from the online Support Center:

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[liofilchem.com/ifu-sds](https://www.liofilchem.com/ifu-sds)

References / Riferimenti

1. US Food and Drug Administration. Bacteriological Analytical Manual (BAM) Chapter 5: Salmonella. Rev 2022-03-18.
2. Public Health England (2022) UK SMI ID 20 issue 4: identification of *Shigella* species. <https://www.gov.uk/government/publications/smi-id-20-identification-of-shigella-species/uk-smi-id-20-issue-4-april-2022-identification-of-shigella-species>.
3. Public Health England (2021) UK SMI ID 24 issue 4: identification of *Salmonella* species. https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/969498/ID_24i4.pdf
4. EN ISO 11133:2014+Amd1:2018+Amd2:2020. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
5. Public Health England (2017) Inoculation of culture media for bacteriology. UK Standards for Microbiology Investigations. Q 5 Issue 2. <https://www.gov.uk/guidance/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
6. BAM Media M61: Hektoen Enteric (HE) Agar. 2017-10-17.
7. ASM (2016) Hektoen Enteric Agar Protocol. American Society for Microbiology.
8. APHA (2015): Compendium of Methods for the Microbiological Examination of Foods. 5th ed. American Public Health Association, Washington, D.C.
9. ISO 21567:2004. Microbiology of food and animal feeding stuffs -- Horizontal method for the detection of *Shigella* spp.
10. Perez J.M., P. Cavalli, C. Roure, R. Renac, Y. Gille, and A. M. Freydiere (2003) Comparison of Four Chromogenic Media and Hektoen Agar for Detection and Presumptive Identification of *Salmonella* Strains in Human Stools. J Clin Microbiol; 41(3):1130–1134.
11. Bisciello N.B. jr. and Schrader J.(1974) Evaluation of Hektoen Enteric Agar for the detection of *Salmonella* in foods and feeds. - Journ of AOAC; 57: 992-996.

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