Wurtz Lactose Agar

Medium for the isolation and differentiation of enterobacteria.

Instructions For Use

ENGLISH

DESCRIPTION

Wurtz Lactose Agar is a solid culture medium used for the isolation and differentiation of enterobacteria from clinical and nonclinical specimens.

This medium allows the differentiation of bacteria on the basis of their ability to utilize lactose as source of energy, resulting in acid production and a color change of the medium to yellow.

TYPICAL FORMULA*

	(g/litre)
Peptone	5.0
Meat Extract	3.0
Sodium Chloride	5.0
Lactose	10.0
Bromothymol Blue	0.075
Agar	15.0

Final pH 7.0 ± 0.2 at 25°C

*Formula may be adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

Peptone and meat extract provide amino acids, nitrogen, carbon, minerals and vitamins for organisms growth. Sodium chloride maintains the osmotic balance of the medium. Lactose is the fermentable carbohydrate. Bromothymol blue is the pH indicator. Agar is the solidifying agent.

While lactose fermenters produce an acidification of the medium, lactose non-fermenting enterobacteria cause an alkalization of the medium as consequence of the decomposition of peptones with production of ammonia.

PREPARATION

<u>Medium in bottles</u>	Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.
--------------------------	---

TEST PROCEDURE

Inoculate the medium by streaking the sample directly over the agar surface to obtain well isolated colonies.

Incubate plates inverted at 35 ± 2°C for 18-24 hours in aerobic atmosphere.

INTERPRETING RESULTS

Lactose fermenters form yellow colonies surrounded by yellow medium.

Non-fermenters grow as blue colonies.

STORAGE

Store bottles and prepared plates at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

Avoid quick temperature shifts of plated medium to prevent condensation.

SHELF LIFE

Medium in bottles: 2 years.

Ready-to-use plates: 6 months.

QUALITY CONTROL

Appearance of Prepared Medium: Slightly opalescent, orange-reddish.

Expected Cultural Response:

Control strain		Inoculum	Incubation	Specification
<i>Escherichia coli</i>	WDCM 00013 (ATCC 25922, NCTC 12241)	≤ 100 CFU	18-24 h/ 35 ± 2°C	Good growth, yellow colonies
<i>Klebsiella pneumoniae</i>	WDCM 00097 (ATCC 13883, NCTC 9633)			
<i>Salmonella</i> Enteritidis	WDCM 00030 (ATCC 13076, NCTC 12694)			Good growth, blue colonies
<i>Shigella flexneri</i>	WDCM 00126 (ATCC 12022, NCTC 12698)			

Please refer to the actual batch related Certificate of Analysis (CoA).

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. For professional use only. Operators must be trained and have certain experience in the laboratory methods. Please read the instructions carefully before using this product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.

Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

The product is available in the configurations listed below. There may be additional product ref. numbers as well. For an updated listing of available products, visit liofilchem.com

Product	Format	Packaging	Ref.
Wurtz Lactose Agar	90 mm Plate	20 plates	10054
Wurtz Lactose Agar	90 mm Plate	100 plates	10054*
Wurtz Lactose Agar	Bottle	6 x 100 ml	442300

This IFU document and the SDS are available from the online Support Center:

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