HELICOBACTER PYLORI AGAR

Selective medium for the isolation of *Helicobacter pylori* from clinical specimens.

TYPICAL FORMULA (g/l)

Pancreatic Digest of Casein	10.0	Vitalex	10.0 ml
Meat Peptic Digest	5.0	Lysed Horse Blood	70.0 ml
Heart Pancreatic Digest	3.0	Amphotericin B	0.005
Yeast Extract	5.0	Cefsulodin	0.005
Corn Starch	1.0	Trimethoprim	0.005
Sodium Chloride	5.0	Vancomycin	0.01
Agar	14.0	Final pH 7.3 ± 0.2	

DESCRIPTION

HELICOBACTER PYLORI AGAR is a selective medium containing Columbia Agar as base and used for the isolation of *Helicobacter pylori* from gastric biopsy specimens.

PRINCIPLE

Peptones provide nitrogen, carbon, amino acids and vitamins. Yeast extract and corn starch are included to supply a carbon energy source in the form of B-complex vitamins. Sodium Chloride maintains the osmotic balance of the medium. Agar is the solidifying agent. Supplementation with vitalex and lysed horse blood provides additional growth factors for fastidious microorganisms. The antimicrobial agents combination inhibits the contaminating flora without reducing the recovery of *H. pylori*.

TECHNIQUE

Inoculate the specimen promptly after collection by streaking it onto the agar surface. Incubate at 35 ± 2°C for 3-5 days in microaerobic atmosphere.

Note: *H. pylori* is extremely sensitive to desiccation and exposure to oxygen. The recovery of the organism from gastric biopsies is improved by direct cultivation as soon as possible after collection. If a delay is expected, transport media such as Gesa Medium (ref. 27001) is recommended.

INTERPRETATION OF RESULTS

Helicobacter pylori cultivates with translucent and non-coalescent colonies, grows at 35°C and does not grow at 25°C, it is urease, catalase and oxidase positive and hippurate negative. Final identification should be done performing appropriate biochemical tests.

STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

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HELICOBACTER PYLORI AGAR

PRESENTATION

Ready-to-use plates (90 mm) containing 22 ± 1 ml of medium

STORAGE

2-8°C

PACKAGING

Ref.	Content	Packaging
10082	20 plates	10 plates in thermally soldered film 2 x 10 plates in cardboard box

pH OF THE MEDIUM

7.3 ± 0.2

USE

HELICOBACTER PYLORI AGAR is a selective medium containing Columbia Agar as base and used for the isolation of *Helicobacter pylori* from gastric biopsy specimens

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Red-brown, opaque medium

SHELF LIFE

75 days

QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
7 days at $22 \pm 1^\circ\text{C}$, in aerobiosis
7 days at $36 \pm 1^\circ\text{C}$, in aerobiosis
- Microbiological control
Inoculum for productivity: 10-100 CFU/ml
Inoculum for selectivity: 10^4 - 10^5 CFU/ml
Incubation Conditions: 24-48 hours at $35 \pm 2^\circ\text{C}$, in microaerobic atmosphere

Microorganism

Helicobacter pylori ATCC® 43504
Escherichia coli ATCC® 25922

Growth

Good
Inhibited

TABLE OF SYMBOLS

LOT	Batch code	IVD	In vitro Diagnostic Medical Device		Manufacturer		Use by		Fragile, handle with care
REF	Catalogue number		Temperature limitation		Contains sufficient for <n> tests		Caution, consult instructions for use		Do not reuse

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LIOFILCHEM® S.r.l.