Mueller Hinton Chocolate Agar

Medium for isolation and cultivation of fastidious bacteria.

INTENDED PURPOSE

Medium for the isolation and cultivation of fastidious organisms, particularly *Neisseria gonorrhoeae* and *Haemophilus influenzae*, from clinical specimens. This medium is intended as an aid in the diagnosis, requiring additional tests to complete the diagnostic results.

DESCRIPTION

Mueller Hinton Chocolate Agar supports satisfactory growth of fastidious strains belonging to the genera *Neisseria* and *Haemophilus*.

TYPICAL FORMULA* (Per Litre of Purified Water)

Tryptone	11.0 g	Cyanocobalamin	0.001 g
Peptone	3.0 g	L-Cystine-HCl	0.02 g
Sodium Chloride	3.0 g	L-Tryptophan	0.02 g
Glucose	2.0 g	Pyridoxine	0.003 g
Starch	1.0 g	D-Pantothenic Acid	0.003 g
Sodium Phosphate Bibasic	2.0 g	NAD	0.003 g
Sodium Acetate	1.0 g	Biotin	0.0003 g
Magnesium Glycerophosphate	0.2 g	Thiamine	0.00004 g
Calcium Gluconate	0.1 g	Adenine	0.01 g
Cobalt Sulphate	0.001 g	Guanine	0.01 g
Copper Sulphate	0.001 g	Xanthine	0.01 g
Zinc Sulphate	0.001 g	Uracil	0.01 g
Ferric Sulphate	0.001 g	Agar	15.0 g
Manganese Chloride	0.002 g	Defibrinated Horse Blood	50.0 ml
Menadione	0.002 g	Vitalex Growth Supplement	10.0 ml
Final pH 7.4 ± 0.2 at 25°C			

*Adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

The medium consists of a nutritive base enriched with several growth factors. The cultivation of *Haemophilus influenzae* and *parainfluenzae* is ensured by NAD which is an essential coenzyme in the respiratory chain of many fastidious microorganisms. *Neisseria gonorrhoeae* requires L-glutamine, which is present in Vitalex growth supplement. Adenine and guanine are used for nucleic acids synthesis of enzymes and coenzymes. The content of p-aminobenzoic acid, which is present in Vitalex growth supplement, is reduced to a minimum in order to reduce the inactivation of sulphonamides and trimethoprim when the medium is used for susceptibility testing of bacteria to these antimicrobics. L-cysteine is a source of organic sulfur and allows the regulation of the redox potential of the medium. Glucose is a source of energy. Horse defibrinated blood is a further growth factor for fastidious microorganisms. Agar is the solidifying agent.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: Test tubes, inoculating loops, swabs, incubator, quality control organisms.

SPECIMENS

Mueller Hinton Chocolate Agar can be used for all types of specimens from infections suspected of containing fastidious organisms, particularly specimens taken from body sites normally sterile (e.g., cerebrospinal fluid, abscesses), and as a subculture medium obtained from blood cultures.

Mueller Hinton Chocolate Agar is mainly used for the non-selective isolation of *Neisseria*, *Haemophilus* and other bacteria that may not grow on routinely used blood agar, such as Columbia agar with 5% sheep blood.

Good laboratory practices for collection, transport and storage of the clinical specimens should be applied.

Refer to specific guidelines for more information about specimen collection and preparation.

TEST PROCEDURE

For isolation of fastidious organisms, streak the specimen directly onto the agar surface immediately after reception in the laboratory. Incubate plates in an inverted position at $36 \pm 1^\circ\text{C}$ for 24-48 hours in an aerobic atmosphere enriched with 5% carbon dioxide (*). It is recommended to inoculate specimens also onto appropriate selective media.

***Note:** Incubation conditions (i.e., temperature, time and atmosphere) may vary depending on the specific requirements of the organism being tested.

INTERPRETING RESULTS

At the end of the incubation period, observe the microbial growth. Identification of the microorganisms isolated must be followed by appropriate additional tests.

STORAGE

Store at $2-8^\circ\text{C}$ away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration. Avoid quick temperature shifts to prevent condensation.

SHELF LIFE

4 months.

QUALITY CONTROL

Appearance of Medium: Opalescent, brown.

Expected Cultural Response:

Control strain		Inoculum	Incubation	Criteria
<i>Haemophilus influenzae</i>	ATCC® 10211	50-100 CFU	18-24 h $36 \pm 1^\circ\text{C}$ 5% CO_2	Good growth
<i>Neisseria gonorrhoeae</i>	ATCC® 43069			
<i>Neisseria gonorrhoeae</i>	ATCC® 49226			
<i>Neisseria meningitidis</i>	ATCC® 13090			

Please refer to the actual batch related Certificate of Analysis (CoA).

PERFORMANCE CHARACTERISTICS

Performance testing of Mueller Hinton Chocolate Agar was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Invalid results can be caused by poor specimen quality, improper sample collection, improper transportation, improper laboratory processing, or a limitation of the testing technology. The operator should understand the principles of the procedures, including its performance limitations, in advance of operation to avoid potential mistakes.

Due to nutritional variation, some strains may result in poor growth or fail to grow on this medium.

In vitro susceptibility of an organism to a specific antimicrobial agent does not necessarily mean that the agent will be effective in vivo. Consult appropriate references for guidance in the interpretation of results.

Mueller Hinton Chocolate Agar is intended as an aid in the diagnosis of infectious diseases, requiring further tests to complete the diagnostic results.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.

- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

See ordering info below. There may be additional product ref. numbers as well. For an updated listing of available products, visit liofilchem.com

Product	Format	Packaging	Ref.
Mueller Hinton Chocolate Agar	Plate 90 mm	20 plates	10335

Revision History

Revision	Release Date	Change Summary
0	2024-03-08	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

This IFU document and the SDS are available from the online Support Center:

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liofilchem.com/ifu-sds

References / Riferimenti

1. Murray PR, Baron EJ, Pfaller MA, Tenover FC, and Tenover FC ed. 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
2. Ericsson HM, and Sherris JC. 1971. Antibiotic sensitivity testing. Report of an international collaborative study. Acta Pathol. Microbiol. Scand. Sec. B, Suppl. 217.
3. Barry AL, Garcia F, and Thrupp LD. 1970. An improved single-disk method for testing the antibiotic susceptibility of rapidly growing pathogens. Am. J. Clin. Pathol. 53: 149-158.
4. Ryan KJ, Schoenknecht FD, and Kirby WMM. 1970. Disc sensitivity testing. Hospital Practice 5:91-100.
5. Bauer AW, Kirby WMM, Sherris JC, and Turck M. 1966. Antibiotic susceptibility testing by a standardized single disk method. Am. J. Clin. Pathol. 45: 493-496.
6. Mueller JH, and Hinton J. 1941. A protein-free medium for primary isolation of the gonococcus and meningococcus. Proc. Soc. Exp. Biol. Med. 48: 330-333.

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