Mueller Hinton II Agar + 2% NaCl

Medium for susceptibility testing of methicillin resistant *Staphylococcus aureus* (MRSA).

TYPICAL FORMULA	(g/l)
Beef Extract	2.0
Acid Hydrolysate of Casein	17.5
Starch	1.5
Sodium Chloride	20.0
Agar	15.0
Final pH 7.3 ± 0.2	

DESCRIPTION

Mueller Hinton II Agar + 2% NaCl is a medium recommended for detection of methicillin resistance in *S. aureus* by using gradient strip method or disk diffusion.

PRINCIPLE

Acid hydrolysate of casein and beef extract supply amino acids, nitrogen, minerals, vitamins, carbon and other nutrients which support the growth of microorganisms. Starch acts as a protective colloid against any toxic molecules which may be present in the medium. Hydrolysis of the starch during autoclaving provides a small amount of glucose that is a source of energy. The presence of a high concentration of NaCl allows the growth of halophile staphylococci. Agar is the solidify agent.

This medium is manufactured to contain low levels of thymine and thymidine and controlled levels of calcium and magnesium ions.

The disk diffusion technique (Kirby-Bauer method) is based on the diffusion, through the agar, of the antimicrobial substance which soaks the paper disc. Each disc has a single concentration of the antimicrobial agent that inhibits the microorganism growth showing an halo around the disc. The diameter of the inhibiting halo is correlated with the Minimal Inhibitory Concentration (MIC).

The gradient strip method is a quantitative assay which uses a predefined concentration gradient of antibiotic applied to the inoculated agar plate for determining the MIC and for resistance detection. After incubation, a symmetrical inhibition ellipse centred along the strip is formed and the MIC, given in µg/mL, is read directly from the scale at the point where the edge of the ellipse intersects the strip.

TECHNIQUE

1. Prepare a standardized suspension of the test organism by the direct colony suspension method.
2. Dip a sterile cotton swab into the adjusted suspension.
3. Inoculate the plate by streaking the swab over the entire agar surface.
4. Apply the disc/strip onto the surface of the inoculated agar plate.
5. Incubate aerobically at 35 ± 2°C for 18-48 hours.

INTERPRETATION OF RESULTS

Read according to the method being used and refer to either CLSI or EUCAST guidelines for interpretation of test results.

STORAGE

10-25°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *in vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. The European Committee on Antimicrobial Susceptibility Testing. Guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance. Version 2.0, 2017.
2. CLSI (2017) Performance Standards for Antimicrobial Susceptibility Testing; 27th ed. CLSI Supplement M100S. Wayne, PA: Clinical and Laboratory Standards Institute.
3. The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. Version 7.0, 2017. <http://www.eucast.org>.
4. The European Committee on Antimicrobial Susceptibility Testing. Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST. Version 7.0, 2017. <http://www.eucast.org>.
5. The European Committee on Antimicrobial Susceptibility Testing. Antimicrobial susceptibility testing - Disk diffusion method. Version 6.0.
6. CLSI M02-A12 (2015) Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard - Twelfth Edition. Wayne, PA: Clinical and Laboratory Standards Institute.
7. Campos, Gill, Hare, Miller (1986) ANTIMICROBIAL AGENTS AND CHEMOTHERAPY. Effect of NaCl Supplementation of Mueller Hinton Broth on Susceptibility of Staphylococci to Aminoglycosides.
8. Bauer et al. (1966) J. Clin. Pathol. 45:493-496.
9. Mueller, J.H., and Hinton (1941) Proc. Soc. Exp. Biol. Med. 48: 330-333.

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Mueller Hinton II Agar + 2% NaCl

PRESENTATION

Ready-to-use plates (90 mm) containing 22 ± 1 ml of medium

STORAGE

10-25°C

PACKAGING

Ref.	Content	Packaging
11206	20 plates	10 plates in thermally soldered film 2 x 10 plates in cardboard box

pH OF THE MEDIUM

7.3 ± 0.2

USE

Mueller Hinton II Agar + 2% NaCl is a medium used for MRSA-testing by either disk or gradient diffusion

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Slightly opalescent, amber

SHELF LIFE

6 months

QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
7 days at $22 \pm 2^\circ\text{C}$, in aerobiosis
7 days at $35 \pm 2^\circ\text{C}$, in aerobiosis
- Microbiological control
Inoculum: McFarland 0.5
Incubation Conditions: 16-20 h at $35 \pm 2^\circ\text{C}$, in aerobiosis

Microorganism		Cefoxitin MTS MIC ($\mu\text{g/mL}$)	Interpretation
<i>Staphylococcus aureus</i>	ATCC® 29213	1-4	Sensible
<i>Staphylococcus aureus</i>	ATCC® 43300	>4	Resistant

TABLE OF SYMBOLS

LOT	Batch code	IVD	In vitro Diagnostic Medical Device		Manufacturer		Use by		Fragile, handle with care
REF	Catalogue number		Temperature limitation		Contains sufficient for >n tests		Caution, consult instructions for use		Do not reuse

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LIOFILCHEM® S.r.l.