MUELLER HINTON AGAR + CLOXACILLIN

Medium dedicated to the confirmation of ESBL-producing Enterobacteriaceae.

TYPICAL FORMULA	(g/l)
Beef Extract	2.0
Acid Hydrolyzed of Casein	17.5
Starch	1.5
Agar	17.0
Cloxacillin	0.25
Final pH 7.3 ± 0.2	

DESCRIPTION

MUELLER HINTON AGAR + CLOXACILLIN is a medium used to increase sensitivity when the confirmation of ESBL-producing Enterobacteriaceae is carried out on isolates with high level expression of AmpC β -lactamases.

ESBLs are enzymes that hydrolyze most penicillins and cephalosporins, inhibited by β -lactamase inhibitors such as clavulanic acid, sulbactam and tazobactam. AmpC-type enzymes hydrolyze penicillins, cephalosporins and monobactam and are poorly inhibited by the classical ESBL inhibitors, especially clavulanic acid.

Phenotypic methods for confirmation of ESBL, such as the combination disc test (CDT), the double-disc synergy test (DDST) test and the ESBL gradient test (Liofilchem® MIC Test Strip ESBL) are based on the inhibition of ESBL activity by clavulanic acid. However, the presence of ESBLs may be masked by an overproduction of AmpC enzyme.

PRINCIPLE

Acid hydrolyzed of casein and meat extract provide amino acids, nitrogen, minerals, vitamins, carbon and other factors which allow the growth of microorganisms. Starch acts as a protective substance against toxic molecules which can be present in the medium. Hydrolysis of starch during sterilization supplies a little amount of glucose which represents a source of energy. Agar is the solidifying agent. Cloxacillin is a good inhibitor of AmpC enzyme.

TECHNIQUE

1. Prepare a standardized suspension of the test organism using either the direct colony suspension or growth method.
2. Dip a sterile cotton swab into the adjusted suspension.
3. Inoculate the surface of the plate by streaking the swab over the entire agar surface.
4. Apply the discs or the strips onto the surface of the inoculated agar plate.
5. Incubate aerobically at 35 ± 2°C for 16-20 hours.

INTERPRETATION OF RESULTS

Combination Disc Test (CDT)

The test is positive if the inhibition zone diameter is ≥ 5 mm larger with clavulanic acid than without.

Double Disc-Synergy Test (DDST)

Positive result is reported when the inhibition zone is increased in the direction the disc containing clavulanic acid.

MIC Test Strip ESBL (MTS-ESBL)

The test is positive if 8-fold reduction is observed in the value of the cephalosporin combined with clavulanic acid compared with the value of the cephalosporin alone or if a phantom zone or deformed ellipse is present

STORAGE

2-8°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

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2-8°C

PACKAGING

Ref.	Content	Packaging
11518	20 plates	10 plates in thermally soldered film 2 x 10 plates in cardboard box

pH OF THE MEDIUM

7.3 ± 0.2

USE

MUELLER HINTON AGAR + CLOXACILLIN is a medium used to increase sensitivity when the confirmation of ESBL-producing Enterobacteriaceae is carried out on isolates with high level expression of AmpC β-lactamases

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Amber medium, slightly opalescent

SHELF LIFE

6 months

QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
7 days at 22 ± 1°C, in aerobiosis
7 days at 36 ± 1°C, in aerobiosis
- Microbiological control
Incubation Conditions: 18-20 hours at 35 ± 2°C

Microorganism		Growth	ESBL phenotype	AmpC phenotype
<i>Escherichia coli</i>	ATCC® 35218	Good	–	–
<i>Klebsiella pneumoniae</i>	ATCC® 700603	Good	+	–
<i>Klebsiella pneumoniae</i>	ATCC® BAA-1144	Good	+	+

TABLE OF SYMBOLS

 LOT	Batch code	 IVD	<i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
 REF	Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult instructions for use	 Do not reuse	

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