Chromatic VRE

Chromogenic medium for the detection of vancomycin resistance in enterococci, directly from clinical specimens.

Instructions For Use

ENGLISH

INTENDED USE

Chromatic VRE is a selective and differential chromogenic medium used for the qualitative and presumptive detection of vancomycin-resistant Enterococci (VRE) directly from clinical specimens.

This medium is intended as an aid in the detection of the following bacteria:

Enterococcus faecalis and *Enterococcus faecium* with acquired vancomycin resistance.

Chromatic VRE is not intended to diagnose infection or guide therapy. Results can be interpreted after incubation for 24 hours. Subculture to a non-selective medium is required for confirming identification, antimicrobial susceptibility testing and epidemiological typing.

SUMMARY AND EXPLANATION

Vancomycin Resistant Enterococci (VRE) were first isolated in Europe in the late 1980s and recently have been associated with nosocomial infections around the world. Acquisition of vancomycin resistance is mediated by several *van* gene clusters. The most commonly encountered vancomycin resistance genes are *vanA* and *vanB*, which are highly transferable between enterococci, mostly *E. faecium* and *E. faecalis*, and can be transferred to other organisms, including *S. aureus* (VRSA).

E. faecium is the most frequently isolated species of VRE in hospitals and typically produces high vancomycin (>128 µg/ml) and teicoplanin (≥16 µg/ml) minimum inhibitory concentrations (MICs). These isolates typically contain *vanA* genes. A *vanB*-containing isolate typically produces lower level resistance to vancomycin (MIC 16 to 64 µg/ml) and is susceptible to teicoplanin (MIC ≤1 µg/ml).

In contrast to *vanA* and *vanB*, *vanC* and *vanD* in VRE are chromosomally associated, i.e. much less transferable. The *vanC* gene generally presents with low-level resistance against vancomycin (MIC 2–32 µg/ml) and teicoplanin (MIC ≤1 µg/ml). Moreover, VRE containing *vanC* are typically *E. gallinarum* or *E. casseliflavus*/*E. flavescens*, and associated with low risk of mortality and rarely cause nosocomial outbreaks. The *vanD* gene, generally, confers moderate-level resistance against vancomycin (MIC 64–128 µg/ml) and teicoplanin (MIC 4–64 µg/ml) and is constitutively expressed.

To limit the spread of glycopeptide resistance between enterococci and to more pathogenic bacteria which in turn can lead to clinical isolates resistant to all antibiotics, several public health institutes now screen for the carriage of *vanA*- and *vanB*-positive VRE in hospitalized patients.

TYPICAL FORMULA*

	(g/litre)
Peptones	30.0
Sodium Chloride	5.0
Chromogenic and Selective Mix	6.5
Agar	15.0

Final pH 7.2 ± 0.2 at 25°C

*Adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

Peptones provide amino acids, nitrogen, minerals, vitamins, and nutrient factors for growth of bacteria. Sodium chloride maintains the osmotic balance of the medium. The chromogenic and selective mix facilitates the identification of bacteria on the basis of the color and colony morphology, inhibiting most of yeasts and bacteria with the exception of VRE. Agar is the solidifying agent.

TEST PROCEDURE

Inoculate the specimen with a swab or loop directly onto the medium, or after enrichment in a broth, and streak for isolation, using a loop. Incubate the agar plates in an inverted position (medium side up) aerobically at 35 ± 2°C for 24 hours. If no colonies are seen, reincubate for a total of 48 hours.

Notes:

Do not inoculate more than one specimen per plate. Chromatic VRE is primarily used with peri-rectal/anal swabs and stool specimens. Avoid exposure to light during incubation as this might destroy the chromogens.

INTERPRETING RESULTS

At the end of incubation, observe the appearance of the colonies and interpret the results as indicated in the ID Table.

ID Table.

Microorganism	Typical colony color
Vancomycin-resistant <i>Enterococcus faecalis</i>	Violet to mauve
Vancomycin-resistant <i>Enterococcus faecium</i>	Blue-green
Other organisms (if not inhibited)	Colorless or natural color

See Appendix I for example pictures.

Notes:

- Some micro-organisms other than enterococci (including yeasts, Gram-negative bacilli, *Pediococcus*) may develop on the medium and produce typical colonies but which generally differ morphologically.
- Growth depends on the requirements of each individual microorganism. A lack of growth or the absence of violet to mauve or blue-green colonies on Chromatic VRE does not preclude the presence of *Enterococcus faecalis* and *Enterococcus faecium* with acquired vancomycin resistance.
- Appropriate susceptibility tests, molecular methods or phenotypic methods must be performed to confirm the presence of VRE isolates.

APPEARANCE

Light amber, slightly opalescent.

STORAGE

Store at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

3 months.

QUALITY CONTROL

To check the performance of the medium the following microbial reference strains can be used.

Microorganism		Inoculum	Incubation	Sepecification
<i>Enterococcus faecalis</i> (vanB)	ATCC® 51299	50-100 CFU	24-48 h / 35 ± 2°C	Good growth, violet to mauve colonies
<i>Enterococcus faecium</i> (vanA)	DSM 13590			Good growth, blue-green colonies
<i>Enterococcus faecalis</i>	ATCC® 29212	10 ⁴ -10 ⁶ CFU		Inhibition
<i>Enterococcus faecium</i>	ATCC® 19434			Inhibition
<i>Escherichia coli</i>	ATCC® 25922			Inhibition

LIMITATIONS

Chromatic VRE is considered for **screening**, and **not for final identification** of VRE.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for professional use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

Product	Format	Packaging	Ref.
Chromatic VRE	Plate 90 mm	20 plates	11621

This document is available from the online Support Center:

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liofilchem.com/ifu-sds

BIBLIOGRAPHY / BIBLIOGRAFIA

- Centers for Disease Control and Prevention. Vancomycin-resistant Enterococci (VRE) in Healthcare Settings. <https://www.cdc.gov/hai/organisms/vre/vre.html>
- Centers for Disease Control and Prevention. Vancomycin-resistant Enterococci (VRE) and the Clinical Laboratory. <https://www.cdc.gov/hai/settings/lab/vreclinical-laboratory.html>
- Flipse J *et al.* (2019) Appearance of vanD-positive Enterococcus faecium in a tertiary hospital in the Netherlands: prevalence of vanC and vanD in hospitalized patients. *Sci Rep.* 9(1):6949.
- Uttley AH, Collins CH, Naidoo J, George RC (1988) Vancomycin-resistant enterococci. *Lancet* 1 (8575-6): 57-8.
- Leclercq R, Derlot E, Duval J, Courvalin P (1988) Plasmid-mediated resistance to vancomycin and teicoplanin in *Enterococcus faecium*. *N. Engl. J. Med.* 319(3): 157-61.
- Chavers LS, Moser SA, Benjamin WH, Banks SE, Steinhauer JR, Smith AM, Johnson CN, Funkhouser E, Chavers LP, Stamm AM, Waiters KB (2003) Vancomycin-resistant enterococci: 15 years and counting. *J. Hosp. Infect.* 53(3): 159-71.
- Courvalin P (2006) Vancomycin Resistance in Gram-positive Cocci. *Clin. Infect. Dis.* 42:S25-34.

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
liofilchem@liofilchem.com

Chromatic VRE


Chromogenic medium for detection of vancomycin-resistant Enterococci (VRE) directly from clinical specimens.

Enterococcus faecium

Enterococcus faecalis

Enterococcus faecium (green-blue)
Enterococcus faecalis (violet to mauve)

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