Chromatic Campylobacter

Chromogenic medium for the presumptive isolation and identification of *Campylobacter* spp.

INTENDED PURPOSE

Chromogenic selective medium for the presumptive isolation and identification of *Campylobacter* directly from clinical and non-clinical specimens. This medium is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Campylobacter is the most frequent cause of gastroenteritis in humans. *C. jejuni*, *C. coli* and *C. lari* are the species most commonly identified as causing human infection. The bacteria live in the intestinal tracts of animals and consequently food of animal origin can become contaminated. Its isolation and identification provide information about the potential health hazard of food.

Chromatic Campylobacter is a chromogenic medium used for the selective isolation and differentiation of *Campylobacter* species from different materials.

Clinical specimens such as stools and rectal swabs can be inoculated directly onto the agar plate.

This medium can also be used for detection of *Campylobacter* in food, animal feed and environmental samples in accordance with the ISO 10272-1.

TYPICAL FORMULA* Per Litre of Purified Water

Peptone	23.2 g
Growth factors	2.01 g
Chromogenic and Selective Mix	1.31 g
Salts	0.6 g
Agar	15.0 g
Final pH 6.9 ± 0.2 at 25°C	

*Adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

Peptone provides amino acids, nitrogen, carbon, minerals and vitamins essential for bacterial growth.

Growth factors provide important nutritive components for bacteria. Salts maintains the osmotic balance of the medium. Chromogenic and Selective Mix allows the identification of bacteria on the basis of the colony color. Agar is the solidifying agent.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: inoculating loops, swabs, incubator, quality control organisms.

SPECIMENS

All types of specimens should be sampled at the acute stage before antimicrobial therapy (where possible) and examined as soon as possible after collection.

Good laboratory practices for collection, transport and storage of the clinical specimens should be applied. Refer to specific guidelines for more information about specimen collection and preparation.

TEST PROCEDURE

Ensure there is no excessive moisture on the plate before use.

Inoculate the plates by directly streaking the specimen onto the agar surface to obtain well-isolated colonies.

Incubate plates in microaerobic conditions at $41.5 \pm 1^{\circ}\text{C}$ for 40-48 hours in an inverted position (agar-side up). Do not extend incubation beyond 48 h as it may increase risk of false positive results.

For more details, consult appropriate guidelines.

INTERPRETING RESULTS

Campylobacter spp. produce pink-red colonies.

Other organisms, if not inhibited, will appear colorless or color different from pink-red.

Notes: Further testing should be conducted to confirm the presumptive identification of *Campylobacter* isolated on this medium (for example Hippurate Hydrolysis). It may be required to subculture to a nonselective nutrient agar before performing appropriate biochemical and serological tests.

STORAGE

Store plates at 2-8°C away from direct light in their original pack until just prior to use, store in the dark. Avoid quick temperature shifts of plated medium to prevent condensation.

Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

3 months

QUALITY CONTROL

Appearance of Prepared Medium: Slightly opalescent, dark amber.

Expected Cultural Response:

Control strain		Inoculum	Incubation	Criteria	Colony color
<i>Campylobacter jejuni</i>	ATCC® 33560	50-100 CFU	40-48 h / 41.5 ± 1°C	Good growth (P _R ≥ 0.5)	Pink-red
<i>Escherichia coli</i>	ATCC® 25922	10 ⁴ -10 ⁶ CFU		Inhibition	---
<i>Enterobacter faecalis</i>	ATCC® 29212				
<i>Candida albicans</i>	ATCC® 10231				

A productivity ratio (P_R) of 0.5 is equivalent to a recovery rate of 50%.

Please refer to the actual batch related Certificate of Analysis (CoA).

PERFORMANCE CHARACTERISTICS

Performance testing of Chromatic Campylobacter was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Invalid results can be caused by poor specimen quality, improper sample collection, improper transportation, improper laboratory processing, or a limitation of the testing technology. The operator should understand the principles of the procedures, including its performance limitations, in advance of operation to avoid potential mistakes.

C. fetus might not grow in this medium.

Disturbing the microaerophilic atmosphere during growth phase may affect recovery of the *Campylobacter* and cause a false negative result. Read plates after the full incubation period for best results.

Chromatic Campylobacter is intended as an aid in the diagnosis of infectious diseases, requiring further tests to complete the diagnostic results. All identification tests should ideally be performed from non-selective agar.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.

- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

See ordering info below. There may be additional product ref. numbers as well. For an updated listing of available products, visit [liofilchem.com](https://www.liofilchem.com)


Product	Format	Packaging	Ref.
Chromatic Campylobacter	Plate 90 mm	20 plates	11644

Revision History

Revision	Release Date	Change Summary
0	2024-03-21	Document creation

This IFU document and the SDS are available from the online Support Center:

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liofilchem.com/ifu-sds

References / Bibliografia

1. Centers for Disease Control and Prevention. 2001. Campylobacter infections. <https://www.cdc.gov/campylobacter/>
2. EN ISO 10272-1:2017. Microbiology of the food chain. Horizontal method for detection and enumeration of Campylobacter spp. Part 1: Detection method
3. EN ISO 11133:2014+Amd1:2018+Amd2:2020. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
4. Epps SV, Harvey RB, Hume ME, Phillips TD, Anderson RC, Nisbet DJ. Foodborne Campylobacter: infections, metabolism, pathogenesis and reservoirs. Int J Environ Res Public Health. 2013 Nov 26;10 12:6292-304.
5. European Food Safety Authority: explains zoonotic diseases Campylobacter.
6. Kaakoush NO, Castaño-Rodríguez N, Mitchell HM, Man SM. Global Epidemiology of Campylobacter Infection. Clin Microbiol Rev. 2015 Jul;28 3:687-720.
7. Murphy H, Cogan T, Humphrey T. Direction of neutrophil movements by Campylobacter-infected intestinal epithelium. Microbes Infect. 2011 Jan;13 1:42-8.
8. Newell DG, Fearnley C. Sources of Campylobacter colonization in broiler chickens. Appl Environ Microbiol. 2003 Aug;69 8:4343-51.
9. Nolan CM, Johnson KE, Coyle MB, Faler K. Campylobacter jejuni enteritis: efficacy of antimicrobial and antimotility drugs. Am J Gastroenterol. 1983 Oct;78 10:621-6.
10. Oyarzabal, O.A., K.S. Macklin, J.M. Barbaree, and R.S. Miller. 2005. Evaluation of Agar Plates for Direct Enumeration of Campylobacter spp. from Poultry Carcass Rinses. Appl. and Environ. Microbio. 71(6):3351-3354.

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