Muller Kauffmann Tetrathionate Novobiocin Broth (MKTn)

Selective enrichment medium for detection of *Salmonella* spp,
according to ISO 6579-1.

DESCRIPTION

Muller Kauffmann Tetrathionate Novobiocin Broth (MKTn) is a liquid medium used for the selective enrichment of *Salmonella* spp, including *S. Typhi* and *S. Paratyphi*.

This medium is recommended by ISO 6579-1 for the examination of food, animal feed and environmental samples from the food production area.

TYPICAL FORMULA

	(g/l)
Meat Extract	4.3
Enzymatic Digest of Casein	8.6
Sodium Chloride	2.6
Calcium Carbonate	38.7
Sodium Thiosulfate	30.5*
Ox Bile	4.78
Brilliant Green	0.0096
Iodine	4.0
Potassium Iodide	5.0
Novobiocin	0.04

Final pH 8.0 ± 0.2 at 25°C

(*) Equivalent to 47.8 g of sodium thiosulfate pentahydrate.

METHOD PRINCIPLE

Meat extract and enzymatic digest of casein provide amino acids, carbon, nitrogen, vitamins and minerals for organisms growth. Sodium chloride maintains the osmotic balance of the medium. Calcium carbonate neutralizes toxic metabolites. Sodium thiosulfate react with iodine and potassium iodide to form tetrathionate, which inhibit the contaminant bacterial flora. Ox bile and brilliant green are additional selective agent used to suppress coliform bacteria and Gram-positive bacteria. Novobiocin is mainly effective against *Proteus* spp.

TEST PROCEDURE

Dilute food sample 1:10 in Buffered Peptone Water, e.g. a 25 g test portion is mixed with 225 ml of BPW. Incubate at 34-38°C for 18 ± 2 h.

Transfer 1 ml of the pre-enrichment culture to a tube containing 10 ml of MKTn broth. Incubate at 37 ± 1°C for 24 ± 3 h. An additional 24 h incubation may be required.

INTERPRETING RESULTS

Turbidity indicates microbial growth.

Use this enrichment culture to inoculate XLD agar plates and a second different selective solid medium. Incubate aerobically at 37 ± 1°C for 24 ± 3 h. Examine for typical colonies of *Salmonella* spp which have a black centre and a lightly transparent zone of reddish color on XLD agar.

STORAGE

Store at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Medium in bottles: 2 years.

Medium in tubes: 1 year.

QUALITY CONTROL

Appearance: Very pale green, opaque, heterogeneous with precipitate.

Expected cultural response:

Control strain	Inoculum	Incubation	Specifications
<i>Salmonella</i> Typhimurium + <i>Escherichia coli</i> + <i>Pseudomonas aeruginosa</i>	WDCM 00031 WDCM 00013 WDCM 00025	24 ± 3 h/ 37 ± 1°C	>10 characteristic colonies on XLD agar
<i>Escherichia coli</i>	WDCM 00013		Partial inhibition, ≤100 colonies on Tryptic Soy Agar
<i>Enterococcus faecalis</i>	WDCM 00009		<10 colonies on TSA

Please refer to the actual batch related Certificate of Analysis (CoA).

WARNING AND PRECAUTIONS

For professional use only. Operators must be trained and have certain experience. Please read the instructions carefully before using this product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.

Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

The product is available in the configurations listed below. There may be additional product ref. numbers as well. For an updated listing of available products, visit [liofilchem.com](https://www.liofilchem.com)

Product	Format	Packaging	Ref.
Muller Kauffmann Tetrathionate Novobiocin Broth (MKTn)	Tube	20 x 10 ml	20072
	Tube	100 x 10 ml	26072
	Plastic Tube	100 x 10 ml	26072P
	Bottle	6 x 100 ml	403110

Revision History

Revision	Release Date	Change Summary
2	2024-10-22	Updated layout and content
1	2018-02-13	Updated ordering info

This IFU document and the SDS are available from the online Support Center:

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