VAGITUBE

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BIBLIOGRAFIA

1. Rayner, C.F.A. 1968. Comparison of culture media for the growth of *Trichomonas vaginalis*. Br. J. Vener. Dis. **44**:63-66
2. Murray, P.R. (ed.) 1999. *Manual of Clinical Microbiology*, 7th ed. American Society of Microbiology, Washington, D.C.
3. Fouts, A.C., Kraus, S.J. 1980. *Trichomonas vaginalis*: reevaluation of its clinical presentation and laboratory diagnosis. J. Infect. Dis. **141**:137-143.

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Prodotto	REF	Σ
VAGITUBE	20340	20 Provette

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F00022

Rev.2 / 12.05.2011

Device for the research by culture of *Trichomonas vaginalis* and yeasts in urogenital samples

DESCRIPTION

VAGITUBE is a device for the search by culture for *Trichomonas vaginalis* and for yeasts present in urogenital samples. The various components of the culture medium include horse serum as a growth supplement and chloramphenicol as a selective supplement to inhibit concomitant bacterial flora. The product can also be used as a medium for maintaining strains of *T.vaginalis* isolated from the samples examined.

KIT CONTENTS

20 Glass test tubes containing 5 mL of liquid medium

- 20 Internal labels
- 1 Instruction sheet

PRINCIPLE OF THE METHOD

VAGITUBE comprises a test tube containing a selective liquid medium, with a formulation that can support the growth of both yeasts and *T. vaginalis*. As selective agent, the medium contains chloramphenicol, a broad-spectrum antibiotic that inhibits both Gram-positive and Gram-negative bacteria.

COMPOSITION

Tryptone	17.0	g/L
Peptone	3.0	g/L
Yeast Extract	10.0	g/L
Maltose	5.0	g/L
L-Cysteine	1.0	g/L
Ascorbic Acid	1.0	g/L
Potassium phosphate monobasic	1.0	g/L
Potassium phosphate dibasic	0.6	g/L
Chloramphenicol	0.1	g/L
Horse serum (inactivated)	100	mL/L
Distilled water	900	mL/L

COLLECTION OF SAMPLES

Samples to be subjected to cultural examination for the search for *T. vaginalis* or yeasts are normally one of the following: vaginal exudate, urinary sediment, prostate secretion or urethral secretion. These materials are collected and then seeded using a swab or pipette.

T. vaginalis is sensitive to changes in temperature. To obtain good sensitivity of the method, it is therefore necessary to inoculate the sample into the culture medium brought either to ambient temperature or better to 36 ± 1 °C, immediately after the sample has been taken and to start incubation immediately.

TEST PROCEDURE

- Bring the number of Vagitubes required to ambient temperature or, better, to 36 ± 1 °C.
- Open the test tubes only immediately before inoculation. Test tubes opened, closed and not used within a few hours must be discarded.
- Inoculate the sample by pressing and rotating the swab against the walls at the bottom of the test tube or, if the sample is liquid, by transferring a few drops with a pipette to the bottom of the test tube.
- Close the test tube again immediately.
- Complete the label with the patient data and date of seeding and glue it to the test tube.
- Incubate at 36 ± 1 °C for 3-7 days.
- Periodically take a drop of medium from the bottom of the test tube with a pipette and examine it under the microscope for the presence of the characteristic motile trophozoites of *T.vaginalis* or of yeast cells.

INTERPRETATION OF THE RESULTS

Normally, there is a growth of *T.vaginalis* and/or yeasts within 3 days, but some authors recommend 7 days of incubation to confirm a sample negative.

T.vaginalis may grow without producing evident signs of turbidity of the culture medium.

Place a drop of medium on a slide, cover it with a cover slide, and examine it immediately under the microscope (100-400 magnification).

T.vaginalis appears in the form of pear-shaped trophozoites, motile thanks to polar flagellae. *Candida spp.* or other yeasts appear as oval or round cells, often with budding.

In the absence of evident signs of microbial growth, a microscope examination is recommended before declaring the sample negative.

The presence of even one trophozoite of *T.vaginalis*, with its characteristic shape and motility, indicates a current infection in the patient.

The presence of yeasts is of pathological significance only in the presence of numerous cells per microscope field after 24 hours of incubation.

QUALITY CONTROL

Each batch of VAGITUBE is subjected to quality control using a culture of *Candida albicans* ATCC 10231 and a culture of *T.vaginalis* ATCC 30001 for positive controls and a culture of *Escherichia coli* ATCC 25922 for the negative control.

PRECAUTIONS

The VAGITUBE kit does not contain substances classified as hazardous under current legislation, nor does it contain harmful substances in concentrations $\geq 1\%$; it therefore does not require a Safety Data Sheet to be available.

VAGITUBE is a disposable device to be used only for *in vitro* diagnostic use. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store VAGITUBE at 2-8 °C in the original packaging. Keep away from sources of heat and avoid excessive changes of temperature. In such conditions the product VAGITUBE will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration such as turbidity of the broth.

DISPOSAL OF USED MATERIAL

After use, VAGITUBE and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

REFERENCES

1. Rayner, C.F.A. 1968. Comparison of culture media for the growth of *Trichomonas vaginalis*. Br. J. Vener. Dis. **44**:63-66
2. Murray, P.R. (ed.) 1999. *Manual of Clinical Microbiology*, 7th ed. American Society of Microbiology, Washington, D.C.
3. Fouts, A.C., Kraus, S.J. 1980. *Trichomonas vaginalis*: re-evaluation of its clinical presentation and laboratory diagnosis. J. Infect. Dis. **141**:137-143.

PRESENTATION

Product	REF	
VAGITUBE	20340	20 Test tubes

TABLE OF SYMBOLS

IVD	In Vitro Diagnostic Medical Device		Do not reuse
REF	Catalogue number		Fragile, handle with care
	Manufacturer		Contains sufficient for <n> tests
	Use by		Caution, consult accompanying documents
	Temperature limitation	LOT	Batch code

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