## Rappaport Vassiliadis Soy (RVS) Broth

Selective liquid medium for detection of *Salmonella* from clinical and nonclinical specimens.

### INTENDED PURPOSE

Liquid medium for selective enrichment of *Salmonella* spp. from clinical specimens, food, animal feed, environmental specimens. This medium conforms to the requirements of ISO 6579-1. This medium is intended as an aid in the diagnosis, requiring additional tests to complete the diagnostic results.

### DESCRIPTION

Rappaport Vassiliadis Soy (RVS) Broth was modified by Vassiliadis by reducing malachite green concentration and increasing incubation temperature, thus offering a better stability of the pH of the prepared medium and optimizing the concentration of magnesium chloride, resulting in an improved recovery of *Salmonellae*. This medium meets the requirements of ISO 6579-1.

### TYPICAL FORMULA\* (Per liter of Purified Water)

|                                                                   |         |
|-------------------------------------------------------------------|---------|
| Enzymatic Digest of Soya                                          | 4.5 g   |
| Sodium Chloride                                                   | 7.2 g   |
| Potassium Dihydrogen Phosphate (KH <sub>2</sub> PO <sub>4</sub> ) | 1.26 g  |
| Dipotassium Hydrogen Phosphate (K <sub>2</sub> HPO <sub>4</sub> ) | 0.18 g  |
| Magnesium Chloride Anhydrous                                      | 13.4 g  |
| Malachite Green                                                   | 0.036 g |
| Final pH 5.2 ± 0.2 at 25°C                                        |         |

\*Adjusted and/or supplemented as required to meet performance specifications.

Grams per liter of purified water

### METHOD PRINCIPLE

Enzymatic digest of soya provides amino acids, nitrogen, carbon, minerals and vitamins for organisms growth. Sodium chloride maintains the osmotic balance of the medium. Potassium phosphates act as a buffer. Magnesium chloride and malachite green are the selective agents. The low pH helps inhibit non-target organisms.

### PREPARATION

#### Dehydrated medium

Suspend 26.6 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Dispense into final containers. Sterilize in autoclave at 115°C for 15 minutes.

### MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: inoculating loops, sterile cotton swabs, physiological solution, incubator, quality control organisms.

### SPECIMENS

Clinical specimens should be sampled at the acute stage, before antimicrobial therapy (where possible) and examined as soon as possible after collection.

Good laboratory practices for collection, transport and storage of the clinical specimens should be applied. Refer to specific guidelines for more information about specimen collection and preparation.

### TEST PROCEDURE

Following the method described in the ISO 6579-1, transfer 0.1 ml of the culture obtained in the pre-enrichment (Buffered Peptone Water, ref. 24099) to a tube containing 10 ml of the RVS Broth.

Incubate at 41.5 ± 1°C for 24 ± 3 hours.

Use the culture obtained to inoculate the selective solid media, i.e. XLD agar (ref. 10056) and a second isolation agar. For more information, see the ISO document.

### Notes:

Before use, allow RVS broth to equilibrate at room temperature if it was stored at a lower temperature. For some products, like dried milk products and cheese, it may be necessary to incubate the selective enrichment medium for an additional 24 h.

After incubation, the selective enrichment can be stored refrigerated at  $5 \pm 3^{\circ}\text{C}$  for up to 72 h.

## INTERPRETING RESULTS

The presence of turbidity compared to uninoculated control, or a pellicle formation indicates microbial growth. Subculture to suitable solid media for complete identification of the isolated colonies.

## STORAGE

The powder is very hygroscopic, store the powder at  $10\text{--}30^{\circ}\text{C}$ , in a dry environment, in its original container tightly closed.

Store tubes at  $10\text{--}25^{\circ}\text{C}$  away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

## SHELF LIFE

Dehydrated medium: 2 years.

Medium in tubes: 2 years.

## QUALITY CONTROL

**Appearance of Dehydrated Medium:** Free-flowing, homogeneous, green.

**Appearance of Prepared Medium:** Clear, blue.

**Expected Cultural Response:**

| Control strains                                                                                                                         | Inoculum              | Incubation               | Specification                                                           |
|-----------------------------------------------------------------------------------------------------------------------------------------|-----------------------|--------------------------|-------------------------------------------------------------------------|
| <i>Salmonella</i> Typhimurium WDCM 00031<br>+ <i>Escherichia coli</i> WDCM 00012 or 00013<br>+ <i>Pseudomonas aeruginosa</i> WDCM 00025 | ≤ 100 CFU             | 24 ± 3 h /<br>41.5 ± 1°C | >10 characteristic colonies<br>on XLD agar or other<br>medium of choice |
| <i>Salmonella</i> Enteritidis WDCM 00030<br>+ <i>Escherichia coli</i> WDCM 00012 or 00013<br>+ <i>Pseudomonas aeruginosa</i> WDCM 00025 |                       |                          |                                                                         |
| <i>Escherichia coli</i> WDCM 00012 or 00013                                                                                             | > 10 <sup>3</sup> CFU |                          | Partial inhibition ≤100<br>colonies on TSA                              |
| <i>Enterococcus faecalis</i> WDCM 00009 or 00087                                                                                        |                       |                          | <10 colonies on TSA                                                     |

Please refer to the actual batch related Certificate of Analysis (CoA).

## PERFORMANCE CHARACTERISTICS

Performance testing of Rappaport Vassiliadis Soy (RVS) Broth was carried out using the QC strains listed above. The results obtained met the established criteria.

## LIMITATIONS

Invalid results can be caused by poor sample quality, improper sample collection, improper transportation, improper laboratory processing, or a limitation of the testing technology. The operator should understand the principles of the procedures, including its performance limitations, in advance of operation to avoid potential mistakes.

Growth depends on the requirements of each individual microorganism. It is therefore possible that certain strains which have specific requirements (substrate, temperature, incubation conditions, etc.) may not develop.

The combination of malachite green, magnesium chloride and low pH does not allow the growth of certain *Salmonella* serovars such as *Salmonella* Typhi and *Salmonella* Paratyphi A.

Colonies of presumptive *Salmonella* must be sub cultured and their identity confirmed by means of appropriate biochemical and serological tests.

## WARNING AND PRECAUTIONS

- 1) For *in vitro* diagnostic use (IVD).
- 2) For laboratory professional use only.

- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

### DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

### BIBLIOGRAPHY

See the references at the end of this document.

### TABLE OF SYMBOLS

See the table of symbols at the end of this document.

### ORDER INFORMATION

| Product                               | Format           | Packaging   | Ref.   |
|---------------------------------------|------------------|-------------|--------|
| Rappaport Vassiliadis Soy (RVS) Broth | Tube             | 20 x 10 ml  | 24400  |
|                                       |                  | 100 x 10 ml | 26400  |
|                                       | Dehydrated media | 500 g       | 610175 |
|                                       |                  | 100 g       | 620175 |

There may be additional product ref. numbers as well. For an updated listing of available products, visit [liofilchem.com](https://www.liofilchem.com)

### Revision History

| Revision | Release Date | Change Summary                                                                           |
|----------|--------------|------------------------------------------------------------------------------------------|
| 0        | 2025-09-15   | Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0 |

In case of malfunctions or defects, contact immediately Liofilchem (\*) or the local representative.

In case of incident associated with the device, notify immediately Liofilchem (\*) or its local representative and the National Competent Authority.

\*Please login to <https://www.liofilchemstore.it/login.php> (user ID and password required) and click on Complaint.

This IFU document and the SDS are available from the online Support Center:

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**[liofilchem.com/ifu-sds](https://www.liofilchem.com/ifu-sds)**

## References / Riferimenti

1. Konforti K, Navon B, Rappaport F. A new enrichment medium for certain salmonellae. J Clin Pathol 1956; 9:261-266.
2. Vassiliadis P, Pateraki E, Papiconomou N, Papadakis J, Trichopoulos D. Nouveau procédé d'enrichissement de salmonella. Ann. Microb. Irist, Pasteur 1976; 127 B: 195.
3. Van Schothorst M, Renaud AM. Dynamics of salmonella isolation with modified Rappaport's medium (R10). J Appl Bacteriol 1983; 54:209
4. ISO 6579-1:2017 Microbiology of the food chain -- Horizontal method for the detection, enumeration and serotyping of Salmonella -- Part 1: Detection of Salmonella spp. - ISO 6579-1:2017/Amd 1:2020 Microbiology of the food chain — Horizontal method for the detection, enumeration and serotyping of Salmonella — Part 1: Detection of Salmonella spp. — Amendment 1: Broader range of incubation temperatures, amendment to the status of Annex D, and correction of the composition of MSRV and SC
5. Baird RM, Corry JEL, Curtis GDW. Pharmacopoeia of Culture Media for Food Microbiology. Proceedings of the 4th International Symposium on Quality Assurance and Quality Control of Microbiological Culture Media, Manchester 4-5 September, 1986. Int J Food Microbiol 1987; 5:254-255.

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