Lowenstein Jensen Medium

Selective medium for the isolation
and cultivation of mycobacteria.

DESCRIPTION

Lowenstein Jensen Medium is a solid medium used for the selective isolation of *Mycobacterium* spp, especially *Mycobacterium tuberculosis*, from clinical specimens as well as for the cultivation of pure cultures of mycobacteria.

TYPICAL FORMULA (per 600 ml purified water)

L-Asparagine	3.6 g
Monopotassium Phosphate	2.5 g
Magnesium Sulfate	0.24 g
Sodium Citrate	0.6 g
Potato Flour	30.0 g
Malachite Green	0.4 g
Glycerol	12 ml
Whole Egg Suspension	1000 ml
Final pH 7.0 ± 0.2 at 25°C	

METHOD PRINCIPLE

Asparagine and potato flour are sources of nitrogen and vitamins. Monopotassium phosphate and magnesium sulfate enhance organisms growth and act as buffers. Glycerol and whole egg suspension provide fatty acids and amino acids required for the metabolism of mycobacteria. The coagulation of the egg albumin during sterilization makes the medium solid. Sodium citrate is converted to citric acid which holds inorganic cations in solution. Malachite green is the selective agent inhibiting the contaminant microbial flora. Malachite green serves as pH indicator as well.

TEST PROCEDURE

Inoculate the medium according to test procedures recommended by the Centers for Disease Control and Prevention (CDC) or consult appropriate references. Protect tubes from light and incubate at 35 ± 2°C for a minimum of 8 weeks in aerobic atmosphere enriched with 5-10% carbon dioxide. To permit the circulation of CO₂ for the initiation of growth, keep caps loosened for the first week of incubation at least. Thereafter, to prevent dehydration, tighten caps; loosen briefly once a week.

INTERPRETING RESULTS

Examine weekly for growth, pigment production and colony morphology. Carry out identification tests according to established laboratory procedures.

N.B. Negative culture results do not rule-out active infection by mycobacteria. Some factors that may lead to unsuccessful cultures are:

- The specimen was not representative of the infectious materials; i.e. saliva instead of sputum.
- The mycobacteria were destroyed during digestion and decontamination of the specimen.
- Gross contamination interfered with the growth of mycobacteria.
- Proper aerobic conditions and increased CO₂ tension were not provided during incubation.

APPEARANCE

Opalescent, viscous, dark blue to green.

STORAGE

Store at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

1 year.

QUALITY CONTROL

Tubes are inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 100-300 CFU.

Inoculum for selectivity: 10^3 - 10^4 CFU.

Incubation conditions: aerobic atmosphere with 5-10%CO₂ at 35 ± 2°C for up to 21 days.

QC Table.

Microorganism		Growth
<i>Mycobacterium intracellulare</i>	ATCC® 13950	Good
<i>Mycobacterium scrofulaceum</i>	ATCC® 19981	Good
<i>Escherichia coli</i>	ATCC® 25922	Partially to completely inhibited

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry, and M.A. Pfaller (2007) Manual of Clinical Microbiology. 9th ed. ASM Press, Washington, D.C.
2. Clinical and Laboratory Standards Institute (2004) Quality Control for Commercially Prepared Microbiological Culture Media; Approved Standard, 3rd ed. M22-A3. CLSI, Wayne, PA
3. Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger, and W.C. Winn, Jr. (1997) Color Atlas and Textbook of Diagnostic Microbiology. 5th ed. Lippincott Williams and Wilkins, Philadelphia, PA.
4. Kent, P.T. and G.P. Kubica (1985) Public Health Mycobacteriology, A Guide for the Level III Laboratory. Dept. of H.H.S. and Centers for Disease Control (CDC), Atlanta, GA.
5. Jensen, K.A. (1932) Zentralbl. Bakteriол. Parasitenkd. Infektionskr. Hyg. Abt. 1 Orig. 125:222-239.
6. Lowenstein, E. (1931) Zentralbl. Bakteriол. Parasitenkd. Infektionskr. Hyg. Abt. 1 Orig. 120:127-129.

PRESENTATION	Format	Packaging	Ref.
Lowenstein Jensen Medium	Flat bottom slant tubes	10 x 9 ml tubes	30118
Lowenstein Jensen Medium	Flat bottom slant tubes	20 x 9 ml tubes	31118
Lowenstein Jensen Medium	Round bottom slant tubes	10 x 8 ml tubes	35000

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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BIBLIOGRAPHIE

1. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry, and M.A. Pfaller (2007) Manual of Clinical Microbiology. 9th ed. ASM Press, Washington, D.C.
2. Clinical and Laboratory Standards Institute (2004) Quality Control for Commercially Prepared Microbiological Culture Media; Approved Standard, 3rd ed. M22-A3. CLSI, Wayne, PA
3. Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger, and W.C. Winn, Jr. (1997) Color Atlas and Textbook of Diagnostic Microbiology. 5th ed. Lippincott Williams and Wilkins, Philadelphia, PA.
4. Kent, P.T. and G.P. Kubica (1985) Public Health Mycobacteriology, A Guide for the Level III Laboratory. Dept. of H.H.S. and Centers for Disease Control (CDC), Atlanta, GA.
5. Jensen, K.A. (1932) Zentralbl. Bakteriol. Parasitenkd. Infektionskr. Hyg. Abt. 1 Orig. 125:222-239.
6. Lowenstein, E. (1931) Zentralbl. Bakteriol. Parasitenkd. Infektionskr. Hyg. Abt. 1 Orig. 120:127-129.

PRESENTATION	Format	Emballage	Ref.
Lowenstein Jensen Medium	Tubes à fond Plat	10 x 9 ml tubes	30118
Lowenstein Jensen Medium	Tubes à fond Plat	20 x 9 ml tubes	31118
Lowenstein Jensen Medium	Tubes à fond Rond	10 x 8 ml tubes	35000

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