LOWENSTEIN JENSEN MEDIUM w/o GLYCEROL

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LOWENSTEIN JENSEN MEDIUM w/o GLYCEROL

Medium for mycobacteria isolation

TYPICAL FORMULA	(g/l)
Magnesium Sulphate	0.24
Asparagine	3.6
Monopotassium Phosphate	2.4
Magnesium Citrate	0.6
Potato Flour	30.0
Malachite Green	0.4
Egg suspension	1000 ml
Final pH 7.0 ± 0.2	

DESCRIPTION

LOWENSTEIN JENSEN MEDIUM w/o GLYCEROL is a medium for mycobacteria isolation from clinical samples or pure cultures.

PRINCIPLE

Eggs suspension provides acids and proteins required for mycobacteria metabolism. The coagulation of egg albumin during sterilization supplies a solid medium for inoculum. Malachite green selectively inhibits contaminants.

TECHNIQUE

1. Inoculate by streaking the sample onto the medium surface.
2. Incubate tubes at 36±1°C for about eight weeks, in atmosphere containing 10% CO₂.
3. Check the growth in the tubes.

INTERPRETATION OF RESULTS

Observe the colonies growth.

STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the storage or transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, do not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for in vitro diagnostic use only and must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. Lowenstein (1931) Zentralbl. Bakteriol. Parasitenkd. Infektionskr. Hyg. Abt. 1 Orig. 120: 127.
2. Jensen (1932) Zentralbl. Bakteriol. Parasitenkd. Infektionskr. Hyg. Abt. 1 Orig. 125: 222.
3. Nolte, F.S., and B. Methcock. (1995) .Mycobacterium, p. 400-437.
4. Musser, J.M. (1995). Antimicrobial resistance in mycobacteria: molecular genetic insights. Clinical Microbiology Reviews 8: 496-514.

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

PRODUCT SPECIFICATIONS

NAME

LOWENSTEIN JENSEN MEDIUM w/o GLYCEROL

PRESENTATION

Slant tubes containing 8.5 ± 0.5 ml of medium

STORAGE

2-8°C

PACKAGING

Ref.	Content	Packaging
30119	10 tubes	10 tubes in cardboard box

pH OF THE MEDIUM

7.0 ± 0.2

USE

LOWENSTEIN JENSEN MEDIUM w/o GLYCEROL is a medium for mycobacteria isolation from clinical samples or pure cultures

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Opaque light green medium

SHELF LIFE

1 year

QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
7 days at $25 \pm 1^\circ\text{C}$, in aerobiosis
7 days at $36 \pm 1^\circ\text{C}$, in aerobiosis
- Microbiological control
Inoculum for productivity: 10-100 UFC/ml
Inoculum for selectivity: 10^4 - 10^5 UFC/ml
Inoculum for specificity: $\leq 10^4$ UFC/ml
Incubation Conditions: up to 8 weeks in atmosphere containing 10% CO_2

Microorganisms

<i>Escherichia coli</i>	ATCC 25922
<i>Mycobacterium fortuitum</i>	ATCC 6841
<i>Mycobacterium intracellulare</i>	ATCC 13950

Growth

Strongly inhibited
Good
Good

TABLE OF SYMBOLS

LOT	Batch code	IVD	In vitro Diagnostic Medical Device		Manufacturer		Use by		Fragile, handle with care
REF	Catalogue number		Temperature limitation		Contains sufficient for <n> tests		Caution, consult instruction for use		Do not reuse

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net