## COLUMBIA AGAR BASE

Medium for fastidious bacteria growth.

### TYPICAL FORMULA (g/l)

|                      |      |
|----------------------|------|
| Peptospecial.....    | 23.0 |
| Starch.....          | 1.0  |
| Sodium Chloride..... | 5.0  |
| Agar.....            | 14.0 |
| Final pH 7.3 ± 0.2   |      |

### DESCRIPTION

COLUMBIA AGAR BASE is a medium which supports the growth of a wide variety of fastidious microorganisms.

### PRINCIPLE

Peptospecial constitutes a source of proteins and free amino acids. Starch is a source of energy for microorganisms. Sodium chloride maintains the osmotic balance of the medium. The addition of SHEEP BLOOD DEFIBRINATED (code 83296) supplies further growth factors for fastidious microorganisms and allows to show hemolytic reactions. Agar is the solidifying agent.

### PREPARATION

Melt the content of one bottle in a boiling water-bath at 100°C (loosing the caps partially unscrewed) until completely dissolved. Cool down to 45-50°C and aseptically add SHEEP BLOOD DEFIBRINATED (code 83296) at the concentration of 5%. Mix well avoiding the formation of bubbles and aseptically distribute into Petri dishes. Allow the medium to solidify. Store the plates in tightly closed containers.

### TECHNIQUE

Spread the inoculum onto the surface of the agar plate and allow to dry. Invert the plates and incubate at 36±1 °C for 24-48 hours.

### INTERPRETATION of RESULTS

Examine plates for growth and hemolytic reactions. Four types of hemolysis on blood agar media can be described:

1. α-hemolysis is the reduction of hemoglobin to methemoglobin in the medium surrounding the colony, causing a greenish discolorization of the medium.
2. β-hemolysis is the lysis of red blood cells, producing a clear zone surrounding the colony.
3. γ-hemolysis indicates no destruction of red blood cells and no change in the color of the medium.
4. δ-hemolysis indicates a partial lysis.

### STORAGE

10-25°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

### WARNING and PRECAUTIONS

The product is not classified as hazardous by current legislation and does not contain harmful substances in concentrations of ≥1%. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

### DISPOSAL of WASTE

Disposal of waste must be carried out according to national and local regulations in force.

### REFERENCES

1. Ellner, P.D., C.J. Stoessel., E. Drakeford, and F. Vasi (1966). A new culture medium for medical bacteriology. *Am. J.Clin. Path.* 45, 502-504.
2. Isenberg, H.D. (ed.) (1992). *Clinical microbiology procedures handbook*, vol. 1 American Society for Microbiology, Washington, DC.

#### LIOFILCHEM Bacteriology Products

64026 ROSETO D.A. (TE) ITALY- Via Scozia- Zona Ind.le

Tel. +39 085 8930745 - Fax +39 085 8930330

Sito Web: <http://www.liofilchem.net> E-Mail: [liofilchem@liofilchem.net](mailto:liofilchem@liofilchem.net)

## PRODUCT SPECIFICATIONS

### NAME

COLUMBIA AGAR BASE

### PRESENTATION

Code 452210

Glass bottles containing 200 ml of medium.

Code 470100

Glass bottles containing 500 ml of medium.

### STORAGE

10-25 °C

### PACKAGE

| Code   | Content            | Packaging        |
|--------|--------------------|------------------|
| 452210 | 6 bottles x 200 ml | 6 bottles in box |
| 470100 | 6 bottles x 500 ml | 6 bottles in box |

### pH OF THE MEDIUM

7.3 ± 0.2

### USE

COLUMBIA AGAR BASE is a medium which supports the growth of a wide variety of fastidious microorganisms.

### TECHNIQUE

Refer to technical sheet of the product.

### APPEARANCE OF THE MEDIUM

Amber medium, opalescent.

### SHELF LIFE

2 years.

### QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control  
7 days at 25 ± 1 °C, in aerobiosis  
7 days at 36 ± 1 °C, in aerobiosis
- Microbiological control  
Inoculum for productivity: 10-100 UFC/ml  
Inoculum for selectivity: 10<sup>4</sup>-10<sup>5</sup> UFC/ml  
Inoculum for specificity: ≤ 10<sup>4</sup> UFC/ml  
Incubation conditions: 48 h at 36 ± 1 °C, in aerobiosis

| Microorganisms                  |            | Growth | Hemolysis |
|---------------------------------|------------|--------|-----------|
| <i>Staphylococcus aureus</i>    | ATCC 25923 | Good   | Beta      |
| <i>Streptococcus pneumoniae</i> | ATCC 6305  | Good   | Alfa      |
| <i>Streptococcus pyogenes</i>   | ATCC 19615 | Good   | Beta      |
| <i>Gardnerella vaginalis</i>    | ATCC 14018 | Good   | Beta      |

### TABLE OF SYMBOLS

| Symbol | Meanings                           |
|--------|------------------------------------|
|        | Catalogue number                   |
|        | In vitro Diagnostic Medical Device |
|        | Manufacturer                       |
|        | Temperature limitation             |
|        | Kit content                        |
|        | Use by                             |
|        | Batch code                         |
|        | Consult accompanying documents     |

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