Urea-arginine screen

Selective test for screening and detection of urogenital mycoplasmas.

INTENDED PURPOSE

Test for screening and detection of urogenital mycoplasmas (*Ureaplasma urealyticum* and *Mycoplasma hominis*) from clinical specimens such as cervico-vaginal swabs, urethral swabs and seminal fluid. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Urea-arginine screen is a test intended for screening and detection of urogenital mycoplasmas directly from clinical specimens.

KIT CONTENT

- 20 Tubes of Urea-arginine screen
- 20 Vials of Physiological solution (7.0 ml)
- 1 Vial of Vaseline oil (3.0 ml)

METHOD PRINCIPLE

The detection of *Mycoplasma hominis* and *Ureaplasma urealyticum* is based on the fermentation of arginine and urea, respectively, which lead to an alkaline pH shift and a broth color change from yellow to red.

Note: The species *Ureaplasma urealyticum* has been divided into two new species: *Ureaplasma parvum* and *Ureaplasma urealyticum*.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: incubator, inoculating loop, pipettes, culture media, quality control organisms.

REAGENTS

- A tube of Urea-arginine screen contains desiccated biochemical substrates and nutritive components, including arginine, urea and a pH indicator
- Physiological solution (0.85% NaCl)
- Vaseline (mineral) oil

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

PREPARATION OF THE CLINICAL SAMPLE

Take a vial of Physiological Solution provided with the kit.

VAGINAL SWAB – URETHRAL SWAB

Immerse the swab (after obtaining the clinical material) in the vial of Physiological Solution and wait 5 minutes. Carefully squeeze the swab against the vial wall so that the clinical material is dispersed homogeneously in the Physiological Solution. Keep the swab dipped in a nutritive broth until the end of the test.

SEMINAL FLUID

Add 0.2 mL of sample to the vial of Physiological Solution, shake and wait 5 minutes before inoculation.

CLINICAL MATERIAL FROM Mycoplasma Transport Broth (Ref. 20158)

Add 1 mL of Mycoplasma Transport Broth, containing the sample under examination, to the vial of Physiological Solution, shake gently and wait 5 minutes before inoculation. Conserve the tube of Mycoplasma Transport Broth until the end of the test.

INOCULATION OF THE SYSTEM

1. Take a tube of Urea-arginine screen and transfer 0.2 mL clinical sample into the tube.
2. Add 2 drops of Vaseline oil to the tube inoculated with the sample.
3. Incubate the tube at $37 \pm 2^{\circ}\text{C}$ for 18-24 hours and up to 48 hours, if needed.

INTERPRETING RESULTS

A positive test is indicated by a broth colour change to red. No color change (the broth remains yellow) indicates that the test is negative.

YELLOW COLOR		MYCOPLASMAS ABSENCE	RED COLOR		MYCOPLASMAS PRESENCE
--------------	---	---------------------	-----------	---	----------------------

STORAGE

2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

2 years.

QUALITY CONTROL

Appearance of Urea-arginine screen: Yellowish powder, yellow solution once reconstituted

Control strain		Incubation	Characteristic reactions
<i>Ureaplasma parvum</i>	ATCC® 27813/7	18-24 h/ $35 \pm 2^{\circ}\text{C}$	Positive reaction: Colour change to red

PERFORMANCE CHARACTERISTICS

Performance testing of Urea-arginine screen was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Test interpretation should not be made prior to 24 hours of incubation. Earlier interpretation may lead to erroneous false-results.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.

- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Urea-arginine screen	20 tests	81079

Revision History

Revision	Release Date	Change Summary
4	2024-09-16	Added notice to report any malfunction, defect or incident.
3	2023-05-15	Updated layout and content in compliance with IVDR 2017/746.

In case of malfunctions or defects, contact immediately Liofilchem (*) or the local representative.

In case of incident associated with the device, notify immediately Liofilchem (*) or its local representative and the National Competent Authority.

*Please login to <https://www.liofilchemstore.it/login.php> (user ID and password required) and click on Complaint.

This IFU document and the SDS are available from the online Support Center:

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