Kovac's Reagent Droppers

Reagent for indole test of Enterobacteriaceae, in dropper bottles.

INTENDED PURPOSE

Reagent for execution of indole test for *Enterobacteriaceae*. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Kovac's Reagent Droppers is used in determining the ability of bacteria, primarily Enterobacteriaceae, to produce indole by the deamination of tryptophan.

KIT CONTENT

- 20 bottles of Kovac's Reagent Droppers (2.5 ml)

METHOD PRINCIPLE

The production of indole has been widely used as an aid in differentiating certain genera and species within the same genus. Bacteria that produce the enzyme tryptophanase are able to degrade the amino acid tryptophan into pyruvic acid, ammonium and indole. The indole produced reacts with p-dimethylaminobenzaldehyde of the Kovac's (indole) reagent to form a red-violet compound. The chief requirement for culturing an organism prior to performing the indole test is that the medium contains a sufficient quantity of tryptophan.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: inoculating loop, pipettes, culture media, physiological solution (0,85 % saline), quality control organisms.

REAGENTS

- 5% (w/v) p-dimethylaminobenzaldehyde dissolved in a solution of 25% hydrochloridric acid and 75% isobutyl alcohol

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Inoculate a pure culture of the organism to be tested into a tube with a tryptophan containing medium, such as peptone water or tryptic soy broth.
2. Incubate tubes at $35 \pm 2^{\circ}\text{C}$ for 24-48 hours.
3. Add 2-3 drops of Kovac's Reagent directly to the tube.

INTERPRETING RESULTS

The formation of a red to purple color ("cherry-red ring") at the surface of the broth within 30 seconds (after addition of Kovac's Reagent Droppers) indicates a positive reaction for indole production. No color change (the reagent layer remains yellow/brown) is considered as a negative result.

STORAGE

2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

2 years.

QUALITY CONTROL

Appearance of Kovac's Reagent Droppers: Yellow-green liquid, clear.

Control strain		Characteristic reactions
<i>Escherichia coli</i>	WDCM 00013 ATCC® 25922 NCTC 12241	Positive reaction: Development of a red ring
<i>Proteus mirabilis</i>	WDCM 00023 ATCC® 29906 NCTC 11938	Negative reaction: Yellow/brown ring

Please refer to the actual batch related Certificate of Analysis (CoA).

PERFORMANCE CHARACTERISTICS

Performance testing of Kovac's Reagent Droppers was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Organisms to be tested must be grown in a tryptophan containing media (for example tryptone water).

Media such as MacConkey Agar should not be used as they contain indicators which could result in the carryover of colour resulting in false positive colour interpretations.

Media containing glucose should not be used as acid production may inhibit indole production. Only pure cultures should be tested as false positive reactions with mixed cultures (of both indole positive and indole negative) are possible.

Allow at least 30 seconds for the color to develop before considering the indole test negative. Pigmented organisms may yield ambiguous results.

Indole produced by certain anaerobes, such as *Clostridium* species, can be rapidly degraded resulting in false negative reactions.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Kovac's Reagent Droppers	20 x 2.5 ml (100-200 tests)	87001

Revision History

Revision	Release Date	Change Summary
1	2024-09-16	Added notice to report any malfunction, defect or incident.
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

In case of malfunctions or defects, contact immediately Liofilchem (*) or the local representative.

In case of incident associated with the device, notify immediately Liofilchem (*) or its local representative and the National Competent Authority.

*Please login to <https://www.liofilchemstore.it/login.php> (user ID and password required) and click on Complaint.

This IFU document and the SDS are available from the online Support Center:

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