Catalase/Oxy Test

Rapid test for the determination of catalase and peroxidase enzymes.

INTENDED PURPOSE

Rapid test for detection of catalase and peroxidase enzymes used for the differentiation and identification of bacterial strain. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Catalase/Oxy Test is a rapid test for the detection of catalase and peroxidase enzymes. Aerobic and facultative anaerobic organisms produce two toxins during normal metabolism, hydrogen peroxide (H₂O₂) and superoxide radical (O₂⁻). These bacteria have two enzymes that detoxify the products of normal metabolism: Catalase and Peroxidase enzymes.

KIT CONTENT

- 30 tubes of Catalase/Oxy Test
- 2 vials of H₂O₂ Reagent (2.5 ml)

METHOD PRINCIPLE

Chlorpromazine is a molecule constituted by an aliphatic chain which represents the substratum for enzymatic activity of catalase and peroxidase. Catalase is capable of converting hydrogen peroxide to water and oxygen.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: inoculating loop, pipettes, culture media, physiological solution (0.85% saline), quality control organisms.

REAGENTS

- A tube of Catalase/Oxy Test contains 0.1 g of desiccated Chlorpromazine.

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Take the number of Catalase/Oxy Test needed from the fridge and allow them to reach room temperature.
2. Add 4 drops of physiological solution (0.85 % saline) to each tube.
3. With a sterile loop pick up one or more than one well isolated colony from a fresh culture plate (18-24 hours old) and emulsify in the tube.
4. Add 3 drops of H₂O₂ Reagent.
5. Observe the almost immediate development of gas bubbles and the foam formation in the tube.

INTERPRETING RESULTS

Observe the turning of the color of the suspension to pink-red in one minute (Peroxidase Test).

Observe the development of bubbles (foam) on the surface of the suspension (Catalase Test).

Identify the microorganism under examination.

STORAGE

At 2-8°C in the original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

1 year.

QUALITY CONTROL

Control strain		Characteristic reactions
<i>Escherichia coli</i>	WDCM 00013 ATCC® 25922 NCTC 12241	Negative peroxidase reaction: Colorless Negative catalase reaction: No formation of bubbles of oxygen
<i>Staphylococcus aureus</i>	WDCM 00034 ATCC® 25923 NCTC 12981	Negative peroxidase reaction: Colorless Positive catalase reaction: Formation of bubbles of oxygen
<i>Listeria monocytogenes</i>	CIP 7836	Negative peroxidase reaction: Colorless Positive catalase reaction: Formation of bubbles of oxygen
<i>Pseudomonas aeruginosa</i>	WDCM 00025 ATCC® 27853 NCTC 12903	Positive peroxidase reaction: color change to pink-red Positive catalase reaction: Formation of bubbles of oxygen

PERFORMANCE CHARACTERISTICS

Performance testing of Catalase/Oxy Test was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

- The most suitable cultures for oxidase testing are those from culture media without dyes, indicators or inhibitors.
- Bacterial colonies taken from media with pH values below 5.5 (eg after carbohydrate metabolism with subsequent acidification of the culture medium) can give a false negative oxidase reaction.
- Colonies taken from media containing nitrates may give unreliable results.
- The use of nickel-chromium loops may cause false positive results due to oxidation of the iron contained in the loops. A platinum or plastic ring or wooden applicator stick is recommended.
- Colonies taken from media containing nitrates may give unreliable results.
- Cytochrome oxidase production may be inhibited by acid generation; False negative reactions can be produced by bacteria grown on media containing fermentable carbohydrates such as MacConkey, TCBS Agar etc.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Catalase/Oxy Test	30 tests	88023

Revision History

Revision	Release Date	Change Summary
1	2024-09-16	Added notice to report any malfunction, defect or incident.
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

In case of malfunctions or defects, contact immediately Liofilchem (*) or the local representative.

In case of incident associated with the device, notify immediately Liofilchem (*) or its local representative and the National Competent Authority.

*Please login to <https://www.liofilchemstore.it/login.php> (user ID and password required) and click on Complaint.

This IFU document and the SDS are available from the online Support Center:

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References / Riferimenti

1. Galeazzi L., Grilli G., Turchetti G., Groppa G., Giunta S.: High Reactivity of the Chlorpromazine/Hydrogen Peroxidase Test with Cytochrome Oxydase-Positive Bacteria. Fems Microbiology Letters 56 (1988) 317-320. Pub. By Elsevier.
2. Murray, Baron, Pfaller, Tenorev and Tenover: Manual of Clinical Microbiology. Current Edition.
3. Edwin, H.Lenette.: Manual of Clinical Microbiology. Current Edition (ASM)

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