Metronidazole Test

Diagnostic discs for presumptive identification of *Gardnerella vaginalis*.

INTENDED PURPOSE

Test used for *Gardnerella vaginalis* identification. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Metronidazole Test are paper discs impregnated with metronidazole and used in a diagnostic workflow for identification of *Gardnerella vaginalis* based on susceptibility to that antimicrobial agent.

KIT CONTENT

- 2 cartridges of 50 discs.

METHOD PRINCIPLE

A characteristic useful for identifying *G. vaginalis* is the presence of a zone of inhibition around a disc containing 50 µg of metronidazole. The mechanism of metronidazole's antibacterial activity is related to the presence of a nitro group in the chemical structure. The nitro group is reduced by a nitroreductase in the bacterial cytoplasm, generating cytotoxic compounds and free radicals that disrupt the host DNA. The selectivity of metronidazole for anaerobic or microaerophilic microorganisms is due to the redox potential of their electron transport components, which are responsible for nitro group reduction and generate toxic metabolites.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: inoculating loop, pipettes, culture media, physiological solution, sterile swabs, 0.5 McFarland turbidity standard, forceps, quality control organisms.

REAGENTS

- Each disc is impregnated with 50 µg of metronidazole.

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Allow product to reach room temperature before use for minimizing condensation on the disc.
2. Prepare a suspension equivalent to a 0.5 McFarland turbidity standard with colonies from a fresh pure culture of the test organism.
3. Using a sterile loop or swab inoculate evenly the surface of a blood agar plate with the adjusted microbial suspension.
4. Using sterile forceps, apply one disc of Metronidazole Test to the surface of the inoculated agar.
5. Incubate plates at $35 \pm 2^\circ\text{C}$ for 18-24 hours under microaerophilic or anaerobic atmosphere.
6. Observe for an inhibition zone around the disc.

INTERPRETING RESULTS

A clear zone of inhibition around the disc indicates that the organism is susceptible (S) to metronidazole.

No zone of inhibition around the disc indicates that the organism is resistant (R).

Gardnerella vaginalis is susceptible to metronidazole, producing an inhibition zone around the disc.

STORAGE

Metronidazole Test must be stored at $-20/+8^\circ\text{C}$ in their original packaging. Do not store them near sources of heat and do not expose them to excessive temperature variations. In such conditions, Metronidazole Test can be used until the expiry date shown on the label. Do not use after this date. Dispose of if they show traces of humidity.

SHELF LIFE

2 years.

QUALITY CONTROL

Control strain		Interpretation
<i>Gardnerella vaginalis</i>	ATCC® 14018	S, zone of inhibition
<i>Enterococcus faecalis</i>	ATCC® 19433	R, no inhibition zone

PERFORMANCE CHARACTERISTICS

Performance testing of Metronidazole Test was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

The diameters of inhibition zones vary between strains and are affected by the inoculum density. Further biochemical and/or serological tests are needed to confirm identification.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Metronidazole Test	2 x 50 discs	9508

Revision History

Revision	Release Date	Change Summary
1	2024-05-13	Added notice to report any malfunction, defect or incident.
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0.

In case of malfunctions or defects, contact immediately Liofilchem (*) or the local representative.

In case of incident associated with the device, notify immediately Liofilchem (*) or its local representative and the National Competent Authority.

*Please login to <https://www.liofilchemstore.it/login.php> (user ID and password required) and click on Complaint.

This IFU document and the SDS are available from the online Support Center: [liofilchem.com/ifu-sds](https://www.liofilchem.com/ifu-sds)

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