

1. Product and Company Identification

Product No.: LAB110

Product Name: Hektoen Enteric Agar

Manufacturer: LAB M Limited
 Topley House
 52 Wash Lane
 Bury
 Lancashire, BL9 6AU
 United Kingdom

T: +44 (0) 161 797 5729
 F: +44 (0) 161 762 9322
 E: labm@idgplc.com
 W: www.idgplc.com

2. Composition and Information on Ingredients

Hazardous ingredients:

Component	CAS number EINECS number	Concentration	Classification	Risk phrases
Sodium Desoxycholate	302-95-4 206-132-7	3.2% w/w	Xi	R 36/37/38
Sucrose	57-50-1 200-334-9	15.8% w/w	-	-
Ammonium Ferric Citrate	1185-57-5 214-686-6	2.0% w/w	-	-

3. Hazard Identification

Most important hazard: Not classified as hazardous.

Health hazard > Inhalation: Moderately irritating to the respiratory system.

Health hazard > Skin: May cause irritation.

Health hazard > Eyes: May cause slight transient irritation.

Health hazard > Ingestion: May cause irritation.

4. First Aid Measures

After inhalation: Remove to fresh air; if breathing becomes difficult seek medical advice.

After skin contact: Wash off with copious amounts of water.

After eye contact: Flush with copious amounts of water; if irritation persists seek medical advice.

After ingestion: Make victim drink plenty of water to dilute stomach contents; seek medical advice.

5. Fire Fighting Measures

Extinguishing media, suitable: Water spray, Carbon dioxide, dry chemical powder or foam.

6. Accidental Release Measures

Personal precautions: Wear disposable gloves, eye protection, dust mask and protective clothing.

Spill: Sweep or vacuum material and place in a designated waste container. Ventilate and use a water rinse to clean the spill site after material pickup is complete.

7. Handling and Storage

Handling: Avoid contact with eyes, skin and clothing. Prevent inhalation of dust.

Storage: Keep container tightly closed. Store in a well ventilated area at 10 to 30°C.

8. Exposure Controls and Personal Protection

Engineering controls: Engineering methods to prevent control or exposure to dust are preferred. Methods include process or personnel enclosure, mechanical ventilation, and control of process conditions.

Protective measures > Hands: Protective disposable gloves.

Protective measures > Eyes: Chemical safety glasses or goggles.

Protective measures > Body: Laboratory coat.

Protective measures > Inhalation: Dust mask.

UK Occupational Exposure Limits (preparation) > *Sucrose*
Long term – 10 mg/M³
Short term – 20 mg/M³

Iron salts (as Fe)
Long term – 1 mg/M³
Short term – 2 mg/M³

9. Physical and Chemical Properties

Physical form: Powder.

Colour: Light green.

pH (at 76 g/L H₂O): 7.3 – 7.7

10. Stability and Reactivity

Stability:	Stable under recommended storage and handling conditions (refer to part 7 of this MSDS). Hygroscopic.
Reactivity:	Avoid strong oxidising agents.

11. Toxicological Information

Acute toxicity:	Quantitative data on the toxicity of this product is not available.
Eye irritation:	The degree of irritation is insufficient to warrant labelling as an irritant.
Skin irritation:	The degree of irritation is insufficient to warrant labelling as an irritant.

12. Ecological Information

Ecotoxic effects:	Quantitative data on the toxicity of this product is not available.
Further ecological data:	No ecological problems are to be expected when the product is handled and used with due care and attention.

13. Disposal Considerations

Product:	Must be disposed of in full compliance with all federal, state and local environmental regulations.
Packaging:	Must be disposed of in full compliance with all federal, state and local environmental regulations.

14. Transport Regulations

UN number:	Not regulated.
ADR/RID classification:	Not regulated for road transport.
IATA classification:	Not regulated for air transport.
IMDG classification:	Not regulated for sea transport.

15. Regulatory Information

Hazard symbol:	Not applicable.
Risk phrases:	Not classified in accordance with EU regulations.
Safety phrases:	Not classified in accordance with EU regulations.



16. Further Information

Reason for Issue: Alteration to format of MSDS.

The information contained in this MSDS is based on the present state of our knowledge. It characterises this product with regards to the appropriate safety precaution. It does not represent a guarantee of the properties of this product. LAB M Limited shall not be held liable for any damage resulting from the handling or from contact with this product.