

**1. Product and Company Identification**

Product No.: LAB148

Product Name: Palcam Agar Base  
(Polymixin, Acriflavine, Lithium Chloride, Ceftazidime, Aesculin, Mannitol Agar Base)

Manufacturer: LAB M Limited  
Topley House  
52 Wash Lane  
Bury  
Lancashire, BL9 6AU  
United Kingdom

T: +44 (0) 161 797 5729  
F: +44 (0) 161 762 9322  
E: [labm@idgplc.com](mailto:labm@idgplc.com)  
W: [www.idgplc.com](http://www.idgplc.com)

**2. Composition and Information on Ingredients**

Hazardous ingredients:

Component	CAS number EINECS number	Concentration	Classification	Risk phrases
Lithium Chloride	7447-41-8 231-212-3	21.1% w/w	Xi	R 36/37

**3. Hazard Identification**

Most important hazard: Irritating to eyes, respiratory system and skin.

Health hazard > Inhalation: Irritating to the respiratory system.

Health hazard > Skin: May cause irritation.

Health hazard > Eyes: May cause slight transient irritation.

Health hazard > Ingestion: Harmful if swallowed.

**4. First Aid Measures**

After inhalation: Remove to fresh air; if breathing becomes difficult, give oxygen. If not breathing, give artificial respiration. Seek medical advice

After skin contact: Wash off with copious amounts of water. If irritation persists, seek medical advice.

After eye contact: Flush with copious amounts of water; if irritation persists seek medical advice.

After ingestion: Make victim drink plenty of water to dilute stomach contents; seek medical advice.

**5. Fire Fighting Measures**

Extinguishing media, suitable: Water spray, Carbon dioxide, dry chemical powder or foam.

---

**6. Accidental Release Measures**

Personal precautions:	Wear disposable gloves, eye protection, dust mask and protective clothing.
Spill:	Sweep or vacuum material and place in a designated waste container. Ventilate and use a water rinse to clean the spill site after material pickup is complete.

---

**7. Handling and Storage**

Handling:	Avoid contact with eyes, skin and clothing. Prevent inhalation of dust.
Storage:	Keep container tightly closed. Store in a well ventilated area at 10 to 30°C.

---

**8. Exposure Controls and Personal Protection**

Engineering controls :	Engineering methods to prevent control or exposure to dust are preferred. Methods include process or personnel enclosure, mechanical ventilation, and control of process conditions.
Protective measures > Hands:	Protective disposable gloves .
Protective measures > Eyes:	Chemical safety glasses or goggles.
Protective measures > Body:	Laboratory coat.
Protective measures > Inhalation:	Dust mask.

---

**9. Physical and Chemical Properties**

Physical form:	Powder.
Colour:	Straw.
pH (at 71 g/L H <sub>2</sub> O):	7.0– 7.4

---

**10. Stability and Reactivity**

Stability:	Stable under recommended storage and handling conditions (refer to part 7 of this MSDS). Hygroscopic.
Reactivity:	Avoid strong oxidising agents .

---

**11. Toxicological Information**

Acute toxicity:	Harmful by ingestion.
Eye irritation:	May cause irritation.
Skin irritation:	May cause irritation.

**12. Ecological Information**

Ecotoxic effects: Quantitative data on the toxicity of this product is not available.

Further ecological data: No ecological problems are to be expected when the product is handled and used with due care and attention.

---

**13. Disposal Considerations**

Product: Hazardous waste. Must be disposed of in full compliance with all federal, state and local environmental regulations.

Packaging: Hazardous waste. Must be disposed of in full compliance with all federal, state and local environmental regulations.

---

**14. Transport Regulations**

UN number: Not classified.

ADR/RID classification: Not regulated for road transport.

IATA classification: Not regulated for air transport.

IMDG classification: Not regulated for sea transport.

---

**15. Regulatory Information**

Hazard symbol: Irritant.

Risk phrases: R 36/37. Irritating to eyes and respiratory system.

Safety phrases: S22. Do not breathe dust.  
S36/37/39. Wear suitable protective clothing, gloves and eye/face protection.

---

**16. Further Information**

Reason for Issue: Alteration to format of MSDS.

---

*The information contained in this MSDS is based on the present state of our knowledge. It characterises this product with regards to the appropriate safety precaution. It does not represent a guarantee of the properties of this product. LAB M Limited shall not be held liable for any damage resulting from the handling or from contact with this product.*