
1. Product and Company Identification

Product No.: LAB014

Product Name: Nutrient Broth No. 2 (B.P.)

Manufacturer: LAB M Limited
Topley House
52 Wash Lane
Bury
Lancashire, BL9 6AU
United Kingdom

T: +44 (0) 161 797 5729
F: +44 (0) 161 762 9322
E: labm@idgplc.com
W: www.idgplc.com

2. Composition and Information on Ingredients

Hazardous ingredients: Dehydrated Culture Media. This product does not contain substances presenting a health hazard as defined by directive 67/548/EEC.

3. Hazard Identification

Most important hazard: Not classified as hazardous.

Health hazard > Inhalation: Moderately irritating to the respiratory system.

Health hazard > Skin: May cause irritation.

Health hazard > Eyes: May cause slight transient irritation.

Health hazard > Ingestion: May cause irritation.

4. First Aid Measures

After inhalation: Remove to fresh air; if breathing becomes difficult seek medical advice.

After skin contact: Wash off with copious amounts of water.

After eye contact: Flush with copious amounts of water; if irritation persists seek medical advice.

After ingestion: Make victim drink plenty of water to dilute stomach contents; seek medical advice.

5. Fire Fighting Measures

Extinguishing media, suitable: Water spray, Carbon dioxide, dry chemical powder or foam.

6. Accidental Release Measures

| | |
|-----------------------|--|
| Personal precautions: | Wear disposable gloves, eye protection, dust mask and protective clothing. |
| Spill: | Sweep or vacuum material and place in a designated waste container. Ventilate and use a water rinse to clean the spill site after material pickup is complete. |

7. Handling and Storage

| | |
|-----------|---|
| Handling: | Avoid contact with eyes, skin and clothing. Prevent inhalation of dust. |
| Storage: | Keep container tightly closed. Store in a well ventilated area at 10 to 30°C. |

8. Exposure Controls and Personal Protection

| | |
|-----------------------------------|--|
| Engineering controls: | Engineering methods to prevent control or exposure to dust are preferred. Methods include process or personnel enclosure, mechanical ventilation, and control of process conditions. |
| Protective measures > Hands: | Protective disposable gloves. |
| Protective measures > Eyes: | Chemical safety glasses or goggles. |
| Protective measures > Body: | Laboratory coat. |
| Protective measures > Inhalation: | Dust mask. |

9. Physical and Chemical Properties

| | |
|----------------------------------|--------------|
| Physical form: | Powder. |
| Colour: | Light straw. |
| pH (at 41 g/L H ₂ O): | 7.1 – 7.5 |

10. Stability and Reactivity

| | |
|-------------|---|
| Stability: | Stable under recommended storage and handling conditions (refer to part 7 of this MSDS). Hygroscopic. |
| Reactivity: | Avoid strong oxidising agents. |

11. Toxicological Information

| | |
|------------------|---|
| Acute toxicity: | Quantitative data on the toxicity of this product is not available. |
| Eye irritation: | May cause irritation. |
| Skin irritation: | May cause irritation. |

12. Ecological Information

| | |
|--------------------------|---|
| Ecotoxic effects: | Quantitative data on the toxicity of this product is not available. |
| Further ecological data: | No ecological problems are to be expected when the product is handled and used with due care and attention. |

13. Disposal Considerations

| | |
|------------|---|
| Product: | Must be disposed of in full compliance with all federal, state and local environmental regulations. |
| Packaging: | Must be disposed of in full compliance with all federal, state and local environmental regulations. |

14. Transport Regulations

| | |
|-------------------------|-----------------------------------|
| UN number: | Not regulated. |
| ADR/RID classification: | Not regulated for road transport. |
| IATA classification: | Not regulated for air transport. |
| IMDG classification: | Not regulated for sea transport. |

15. Regulatory Information

| | |
|-----------------|---|
| Hazard symbol: | Not applicable. |
| Risk phrases: | Not classified in accordance with EU regulations. |
| Safety phrases: | Not classified in accordance with EU regulations. |

16. Further Information

| | |
|-------------------|-------------------------------|
| Reason for Issue: | Alteration to format of MSDS. |
|-------------------|-------------------------------|

The information contained in this MSDS is based on the present state of our knowledge. It characterises this product with regards to the appropriate safety precaution. It does not represent a guarantee of the properties of this product. LAB M Limited shall not be held liable for any damage resulting from the handling or from contact with this product.