

Safety Data Sheet LAB022 Reinforced Clostridial Medium Issue: 17/Dec/2004

6. Accidental Release Measures

Personal precautions: Wear disposable gloves, eye protection, dust mask and protective clothing.

Spill: Sweep or vacuum material and place in a designated waste container. Ventilate and use a water rinse to clean the spill site after material pickup is complete.

7. Handling and Storage

Handling: Avoid contact with eyes, skin and clothing. Prevent inhalation of dust.

Storage: Keep container tightly closed. Store in a well ventilated area at 10 to 30°C.

8. Exposure Controls and Personal Protection

Engineering controls: Engineering methods to prevent control or exposure to dust are preferred. Methods include process or personnel enclosure, mechanical ventilation, and control of process conditions.

Protective measures > Hands: Protective disposable gloves.

Protective measures > Eyes: Chemical safety glasses or goggles.

Protective measures > Body: Laboratory coat.

Protective measures > Inhalation: Dust mask.

9. Physical and Chemical Properties

Physical form: Powder.

Colour: Light straw.

pH (at 38 g/L H₂O): 6.6 – 7.0

10. Stability and Reactivity

Stability: Stable under recommended storage and handling conditions (refer to part 7 of this MSDS). Hygroscopic.

Reactivity: Avoid strong oxidising agents.

11. Toxicological Information

Acute toxicity: Quantitative data on the toxicity of this product is not available.

Eye irritation: May cause irritation.

Skin irritation: May cause irritation.
