


# Liofilchem®

# CONTACT SLIDE

flexible contact slide for microbiological screening of surfaces

ISO 18593



Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.

printed in Italy

© Liofilchem 2025

Rev. 9 - 3.09.2025

ref. 6550616

Description	Packaging	# Product with minimum order Ref.
<b>Contact Slide Chrom 1</b> (Chromatic™ E.coli-coliform / Chromatic™ Salmonella) Chromogenic selective media for detecting <i>E. coli</i> , coliforms and <i>Salmonella</i> spp. on surfaces or in food and waters.	20 slides	525292
<b>Contact Slide Chrom 2</b> (Chromatic™ E.coli-coliform / PCA+TTC + Neutralizing) Chromogenic selective medium for detecting <i>E. coli</i> , coliforms and for total bacterial count.	20 slides	525392
<b>Contact Slide Chrom 3</b> (Chromatic™ E.coli-coliform / Baird Parker Agar + Neutralizing) Chromogenic selective medium for detecting <i>E. coli</i> , coliforms and <i>Staphylococcus aureus</i> with inactivation of disinfectants.	20 slides	525462
<b>Contact Slide 1</b> (PCA + TTC + Neutralizing / V.R.B.G. Agar + Neutralizing) Total bacterial count, enterobacteria, with inactivation of disinfectants.	20 slides 120 slides	525262 53526
<b>Contact Slide 2</b> (PCA + TTC + Neutralizing / Rose Bengal CAF Agar + Neutralizing) Total bacterial count, yeasts and moulds, with inactivation of disinfectants.	20 slides 120 slides	525272 53527
<b>Contact Slide 3</b> (PCA + TTC + Neutralizing / Vogel-Johnson Agar) Total bacterial count with inactivation of disinfectants, staphylococci.	20 slides 120 slides	525302 53530
<b>Contact Slide 4</b> (Cetrimide Agar + Neutralizing / Rose Bengal CAF Agar + Neutralizing) <i>Pseudomonas</i> spp., yeasts and moulds.	20 slides 120 slides	525282 53528
<b>Contact Slide 5</b> (VRBG Agar + Neutralizing / Bile Esculin Agar + Neutralizing) Enterobacteria, fecal streptococci, with inactivation of disinfectants.	20 slides 120 slides	525312 53531
<b>Contact Slide 8</b> (PCA + TTC + Neutralizing / V.R.B.L. Agar + Neutralizing) Total bacterial count, coliforms, with inactivation of disinfectants.	20 slides 120 slides	525622 53562
<b>Contact Slide 9</b> (PCA + TTC + Neutralizing / PCA + TTC + Neutralizing) Total bacterial count with inactivation of disinfectants.	20 slides 120 slides	525322 53532
<b>Contact Slide 10</b> (Rose Bengal CAF Agar + Neutralizing / Vogel-Johnson Agar) Yeasts and moulds with inactivation of disinfectants, staphylococci.	20 slides 120 slides	525452 # 53545 #
<b>Contact Slide TSA + Neutralizing / TSA + Neutralizing</b> Total bacterial count with inactivation of disinfectants.	20 slides	525472
<b>Contact Slide TSA + Neutralizing / Rose Bengal CAF Agar + Neutralizing</b> Total bacterial count with inactivation of disinfectants, yeasts and moulds.	20 slides 120 slides	525482 53548
<b>Contact Slide PCA + TTC + Neutralizing / Dichloran Rose Bengal CAF + Neutralizing</b> Total bacterial count with inactivation of disinfectants, yeasts and moulds.	120 slides	53543
<b>Contact Slide VRBG Agar / TSA</b> Enterobacteria, Total bacterial count.	20 slides	525382



**Liofilchem srl**  
Headquarters, Manufacturing Site,  
International distribution  
Roseto degli Abruzzi, Italy

**Liofilchem, Inc.**  
US Distribution Center  
Waltham, MA, USA

**Liofilchem Nordics ApS**  
Scandinavia Customer Center  
Copenhagen, Denmark

