



HALALSELECT™

Sampling Devices



Halal

All Halal Select™ Sampling Devices have received halal certification from the American Halal Foundation. They are available pre-hydrated with our proprietary HiCap™ Neutralizing Broth, which contains no animal derived ingredients, and is designed specifically for use in food production environments.



Halal Select™ Sponge Sampler

Ideal for collecting environmental surface samples without touching the sponge. The innovative twist-off handle holds the sponge securely but is quick and easy to remove after collection.



US Patent 8,127,627 B2

Item #	Description
HS-SPONGE-10HC	Halal Select™ Polyurethane Sponge Sampler with 10 ml HiCap™ Neutralizing Broth (100 per case)
HS-SPONGE-DRY	Halal Select™ Polyurethane Sponge Sampler, dry (100 per case)



Halal Select™ Swab Sampler

The polyurethane foam tip and strong, flexible shaft, combined with the easy-to-grip large cap, allows users to press firmly and vigorously scrub surfaces

Item #	Description
HS-SWAB-HC-P	Halal Select™ Swab Sampler, Large Tip with 1 ml HiCap™ Neutralizing Broth (100 per case)
HS-SWAB-5HC	Halal Select™ Swab Sampler, Large Tip with 5 ml HiCap™ Neutralizing Broth (100 per case)
HS-SWAB-10HC	Halal Select™ Swab Sampler, Large Tip with 10 ml HiCap™ Neutralizing Broth (100 per case)
HS-SWAB-DRY	Halal Select™ Swab Sampler, Large Tip, Dry (100 per case)
HS-SWAB-1HC-S	Halal Select™ Swab Sampler, Small Tip with 1 ml HiCap™ Neutralizing Broth (100 per case)



US Patent 9,027,420



www.worldbioproducts.com

World Bioproducts, LLC
Phone: 877-260-6441 International: +1 425-242-4153
Email: intlcare@worldbioproducts.net