



Reference A04137

One pair SteriSox® jersey, 15 ml tryptone salt in twirl type sampling bag, bulk packed

Standard ISO 13307
Microbiology of food and animal feed
Primary production stage
Sampling techniques



European Union regulations

Use:

Sterile surface sampling kit in poultry farming for Salmonella monitoring. For collecting dust and droppings on the ground, in the production environment of species concerned by the European regulations: *Gallus gallus* and turkeys. Suggested associated product: reference A04050, sterile package of one pair of overshoes and one pair of gloves.

Features:

- **SteriSox® jersey**
 - 1 pair of shoes made of knitted jersey
 - With elastic band
- **Diluent**
 - Tryptone salt
 - Calibrated dose : 15ml
- **Sampling bag**
 - Easy to open, twirl type, 18 x 32 cm
 - Writing area for the sample information



Reference	Custom code	Packaging	Cardboard box size
A04137	6307901000	200 units per box (8 zip bag packaging, size 30 x 40 cm of 25 units)	60 x 40 x 25 cm

Sterility and traceability:

Sterilization	Sterilization certificate	Shelf life	Traceability *
Sterilization by gamma rays Calibrated dose minimum 20 kiloGrays	in each box	12 months	On each product/box



* Printed information: reference, designation, batch number, diluent code, production and expiry date