



SteriCloth

Reference A04601

**Sterile swab cloth, 32 x 17 cm, dry
in blue sampling bag with laser pre-cut**



ISO 18593 standard
horizontal methods for sampling
techniques on surfaces

Use:

Sterile sampling device adapted to sample work surfaces in food industries.
Pathogens detected: Listeria, Salmonella, other ones.

Features:

- **Swab cloth**
 - Blue, unwoven, honeycombed structure
 - 32 x 17 cm
- **Diluent**
 - Absence (dry)
- **Blue sampling bag**
 - Easy to open, twirl type, 18 x 30 cm
 - Large writing area for the sample information



Reference	Custom code	Packaging	Cardboard box size
A04601	5603139090	200 units per box (10 zip bag packaging, size 30 x 40 cm of 20 units)	40 x 30 x 30 cm

Sterility and traceability:

Sterilization	Sterilization certificate	Shelf life	Traceability *
Sterilization by gamma rays Calibrated dose minimum 20 kiloGrays	in each box	24 months	On each product/box



* Printed information: reference, designation, batch number, diluent code, production and expiry date

31/10/2023

