



SteriCloth

## Reference A04613

**Sterile swab cloth, 34 x 37 cm, tryptone salt at 10% neutralizing agent, in blue sampling bag with laser pre-cut**



**ISO 18593 standard**  
horizontal methods for sampling techniques on surfaces

### Use:

Sterile sampling device adapted to sample work surfaces in food industries. Used to check the quality of the disinfection processes, even in presence of residual disinfectants which may be present on the sampled surface. Pathogens detected: Listeria, Salmonella, other ones.

### Features:

- **Swab cloth**
  - Blue, unwoven, honeycombed structure
  - 34 x 37 cm
- **Diluent**
  - Tryptone salt at 10% neutralizing agent
  - Calibrated dose: 25 ml
- **Blue sampling bag**
  - Easy to open, twirl type, 18 x 30 cm
  - Large writing area for the sample information



Reference	Custom code	Packaging	Cardboard box size
A04613	5603139090	140 units per box (7 zip bag packaging, size 30 x 40 cm of 20 units)	40 x 30 x 30 cm

### Sterility and traceability:

Sterilization	Sterilization certificate	Shelf life	Traceability *
Sterilization by gamma rays Calibrated dose minimum 20 kiloGrays	in each box	12 months	On each product/box



Before



After

\* Printed information: reference, designation, batch number, diluent code, production and expiry date

31/10/2023