



SteriSox® knitted jersey

To control *Salmonella* in poultry farms

• Sterility

Gamma-ray sterilization by ionization at 20 Kilograys minimum. Sterility of each batch of product verified by an external COFRAC-certified laboratory.

• Boot swab quality

Ensures efficient collection of droppings and dust in livestock buildings. Once in the laboratory, **SteriSox® knitted jersey** with excellent absorption capacity and hydrophilic properties ensure optimum release of germs in the pre-enrichment phase.

• Diluents

Our diluents are manufactured to the ISO 11133 standard.

• Bag quality

Our sampling bags are ideally suited for lab blending, and feature easy opening and resealing to ensure sample integrity throughout transport.

Kit contents:

• Sampling media

jersey boot swab with elasticated top that can be worn over any type of boots, including large-sized ones.

• The diluent and neutralizing agent

are formulated pursuant to standards; constant, calibrated volume for all kits.

• A sampling bag

with integrated resealable zip, and a write-on area to identify the sample. Seals tested and approved for use in a mixer.

• Sterile gloves (optional)

ambidextrous, latex-free.

• Over-shoes (optional)

strong plastic with elasticated top.



European
regulation



ISO 13307 standard

Product guarantee

- Validity guaranteed up to 1 year for pre-moistened products, 2 years for dry products
- Products are delivered with the certificate of ionization
- Traceability by batch number specified on each kit, box and certificate of ionization

Options available:

- Number of boot swab
- Diluent
- Pair of gloves
- Over-shoes





Ref. A04136 and A04137



Ref. A04139



Ref. A04142 and A04148



Ref. A04050S

During breeding

Reference	Pair of SteriSox®	Diluent	Volume (in ml)	Bag marking colour	Number of bag	Pair of gloves +2 over-shoes	Pack	Quantity per box
A04130	1	Dry	-	●	1	No	Of 25	600
A04131	1	Dry	-	●	1	Yes	Unitary	200
A04136	1	Tryptone salt	15	●	1	Yes	Unitary	75
A04137	1	Tryptone salt	15	●	1	No	Of 25	200
A04142	2	Tryptone salt	30	●	1	No	Of 15	120
A04148	2	Tryptone salt	30	●	1	Yes	Unitary	50
A04145	5	Tryptone salt	70	●	2	Yes	Unitary	80
A04149	5	Tryptone salt	70	●●	2	No	Unitary	50

After cleaning

Reference	Pair of SteriSox®	Diluent	Volume (in ml)	Bag marking colour	Number of bag	Pair of gloves +2 over-shoes	Pack	Quantity per box
A04139	1	Peptone water at 10% neutralizing agent	15	●	1	No	Of 25	50

Accessory kit 2 gloves + 2 over-shoes

Reference	Glove colour	Over-shoe colour	Quantity per box
A04050S	Blue	Transparent	100

Scan the QR Code to discover all the SteriSox®



Contact our sales representative

✉ commercial@sodibox.com ☎ +33 (0)2 98 06 62 20

Discover all our products on

www.sodibox.com

