

Technical Specification Sheet



Meat Peptone No. 3 SKU: 700003597, 700004702, 700003598, 700004703 NCM0246

Intended Use

Meat Peptone No. 3 is for use in preparing microbiological culture media in a laboratory setting, and is not intended for use in the diagnosis of disease or other conditions in humans.

Description

Meat Peptone No. 3 is a highly nutritious enzymatic digest of meat and may be used as an ingredient in culture media.

Precaution

Refer to SDS

Test Procedure

Refer to appropriate references for specific procedures using Meat Peptone.

Quality Control Specifications

PHYSICAL/CHEMISTRY:	EXPECTED RESULTS
Appearance of Powder:	Light beige
Solubility in water at 1%	Total
Clarity	Clear, yellow color
pH of 1% solution	6.7 - 7.3
Total Nitrogen	11.0 - 13.0%
Amino Nitrogen	2.5 - 5.0%

Expiration

Refer to expiration date stamped on container. Meat Peptone No. 3 should be discarded if not free flowing, or if the appearance has changed from the original color. Expiry applies to Meat Peptone in its intact container when stored as directed.

Storage

Store powder at 2-30°C away from direct sunlight. Once opened and recapped, place container in a low humidity environment at the same storage temperature. Protect from moisture and light by keeping container tightly closed.

Record ID: FS-TSS-0378 Revision Number: 3.0 Effective Date: 2023-11-17 12:00 AM EST



620 Leshar Place • Lansing, MI 48912
800-234-5333 (USA/Canada) • 517-372-9200
foodsafety@neogen.com • foodsafety.neogen.com