

**Terrific Broth, Modified**  
**SKU: 700003640, 700003641, 700003642, 700003643**  
**NCM0278**

### Intended Use

Terrific Broth, Modified is used with glycerol in cultivating recombinant strains of *Escherichia coli* in a laboratory setting. Terrific Broth, Modified is not intended for use in the diagnosis of disease or other conditions in humans.

### Description

Terrific Broth, Modified was developed by Tartoff and Hobbs to increase yields of plasmid DNA from transformed *E. coli* and is a popular choice of media available for plasmid preps. While cells grown in Terrific Broth, Modified can grow slower than other liquid media, the final cell density and plasmid yields are 2 – 4 fold higher. The higher yield is a result of the increased concentration of Enzymatic Digest of Casein and Yeast Extract in this formula. Glycerol is used as a carbohydrate source, and unlike glucose, is not fermented to acetic acid.

### Typical Formulation

Yeast Extract	24.0 g/L
Enzymatic Digest of Casein	12.0 g/L
Dipotassium Phosphate	9.4 g/L
Monopotassium Phosphate	2.2 g/L
Final pH: 7.2 ± 0.2 at 25°C	

### Supplement

Glycerol, 8 mL

Formula is adjusted and/or supplemented as required to meet performance specifications.

### Precaution

Refer to SDS

### Preparation

1. Dissolve 47.6 g of the medium in one liter of purified water containing 8 mL of glycerol.
2. Mix thoroughly.
3. Autoclave at 121°C for 15 minutes.

### Quality Control Specifications

**Dehydrated Appearance:** Powder is homogeneous, free-flowing, and light beige.

**Prepared Appearance:** Prepared medium is clear to trace hazy and light to medium amber, with none to light precipitate.

**Expected Cultural Response:** Cultural response in Terrific Broth, Modified at 35 ± 2°C and examined for growth after 18 - 24 hours incubation.

# Technical Specification Sheet



Record ID: FS-TSS-0337 Revision Number: 1.0 Effective Date: 2023-08-29 12:00 AM EDT

Microorganism	Approx. Inoculum (CFU)	Expected Results
<i>Bacillus subtilis</i> ATCC® 6051	10 - 300	Growth
<i>Escherichia coli</i> ATCC® 25922	10 - 300	Growth
<i>Escherichia coli</i> ATCC® 23724	10 - 300	Growth
<i>Escherichia coli</i> ATCC® 33694	10 - 300	Growth
<i>Escherichia coli</i> ATCC® 33849	10 - 300	Growth
<i>Escherichia coli</i> ATCC® 39403	10 - 300	Growth
<i>Escherichia coli</i> ATCC® 47014	10 - 300	Growth
<i>Escherichia coli</i> ATCC® 53868	10 - 300	Growth

The organisms listed are the minimum that should be used for quality control testing.

## **Test Procedure**

Consult appropriate references for recommended test procedures.

## **Results**

After sufficient incubation, growth is exhibited as turbidity in the medium.

## **Expiration**

Refer to expiration date stamped on the container. The dehydrated medium should be discarded if not free flowing, or if appearance has changed from the original color. Expiry applies to medium in its intact container when stored as directed.

## **Limitation of the Procedure**

Due to nutritional variation, some strains may be encountered that grow poorly or fail to grow on this medium.

## **Storage**

Store dehydrated culture media at 2-30°C away from direct sunlight. Once opened and recapped, place container in a low humidity environment at the same storage temperature. Protect from moisture and light by keeping container tightly closed.

## **References**

1. Tartoff, K. D., and C. A. Hobbs. 1987. Improved media for growing plasmid and cosmid clones. Bethesda Research Laboratories Focus 9:12.
2. Sambrook, J., E. F. Fritsch, and T. Maniatis. 1989. Molecular cloning: a laboratory manual, 2<sup>nd</sup> ed. Cold Spring Harbor Laboratory, Cold Spring Harbor, N. Y.

