

Lactic Acid 10% (NCM4011)

Intended Use

Lactic Acid 10% is for the isolation of Gram Positive bacteria. Lactic Acid 10% is intended for use in culture media for the microbiological analysis of food, animal feed, water and environmental samples, and is for use with:

| | |
|---------|----------------------|
| NCM0018 | Potato Dextrose Agar |
| NCM0093 | Malt Extract Agar |

Typical Formulation by Vial Size

| | Per 1 mL Vial NCM4011-1 |
|-------------|----------------------------|
| Lactic Acid | 100 mg |

Precaution

Refer to SDS

Preparation

NCM4011 is ready to use.

Quality Control Specifications

Appearance: Colorless, clear solution.

Storage

Store sealed vial at 2-8°C away from direct sunlight.

Expiration

Refer to expiration date stamped on the container. Expiry applies to medium in its intact container when stored as directed.

Precautions

Discard if supplement demonstrates any precipitate.